

**Method Comparison Study Report for the ISO 16140-2:2016 validation of
Brilliance™ Staph 24 Agar for the enumeration of coagulase positive
Staphylococci species in a broad range of foods**

MicroVal study number: 2008LR11 renewal

Method/Kit name: Brilliance™ Staph 24 Agar

Report version: MCS ILS summary report 28/03/2019

MicroVal Expert Laboratory: Campden BRI (Linda Everis and Gail Betts gail.betts@campdenbri.co.uk)

Foreword

This report is prepared in accordance with ISO 16140-2:2016 and MicroVal technical committee interpretation of ISO 16140-2 v.1.0

Company: Thermo Fisher Scientific
Wade Road,
Basingstoke,
Hampshire,
UK

Expert Laboratory: Campden BRI

Method/Kit name: Brilliance™ Staph 24 Agar

Validation standard: ISO 16140-2:2016 Microbiology of the food chain —Method validation —Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method

Reference methods:

ISO 6888-1:1999 Incorporating Amendment No 1 and corrigendum No 1. Microbiology of food and animal feeding stuffs- Horizontal method for the enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species) —Part 1: Technique using Baird-Parker agar medium

ISO 6888-1:1999 DAM2:2017(E). Microbiology of food and animal feeding stuffs- Horizontal method for the enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species) —Part 1: Technique using Baird-Parker agar medium. AMENDMENT 2 Inclusion of an alternative confirmation procedure

Scope of validation: A broad range of foods based on categories

1. Dairy products
2. Fishery products,
3. Chocolate, bakery and confectionery products,
4. Meat and meat products
5. Multi-component foods

Certification organisation: Lloyd's Register

List of abbreviations

- AL	Acceptability Limit
- AP	Accuracy Profile
- Art. Cont.	Artificial contamination
- CFU	Colony Forming Units
- CL	confidence limit (usually 95%)
- EL	Expert Laboratory
- \bar{D}	Average difference
- g	Gram
- h	Hour
- ILS	Interlaboratory Study
- Inc/Ex	Inclusivity and Exclusivity
- LOQ	Level of Quantification
- MCS	Method Comparison Study
- min	minute
- ml	Millilitre
- MR	(MicroVal) Method Reviewer
- MVTC	MicroVal Technical Committee
- EL	Expert Laboratory
- n	number of samples
- na	not applicable
- neg	negative (target not detected)
- NG	no growth
- nt	not tested
- RT	Relative Trueness
- SD	standard deviation of differences
- 10 ⁻¹ dilution	10-fold dilution of original food
- 10 ⁻² dilution	100-fold dilution of original food
- PSD	Peptone salt diluent
- BSA	Brilliance™ Staph 24 Agar
- BPA	Baird parker Agar
- RPFA	Rabbit plasma fibrinogen Agar
- BHIB	Brain Heart Infusion Broth

Contents

1	Introduction	5
2	Method protocols	6
2.1	Reference method	6
2.2	Alternative method	6
2.3	Study design	7
3	Method comparison study	7
3.1	Relative trueness study	7
3.1.1	Number of samples	7
3.1.2	Test sample preparation	9
3.1.3	Protocols applied during the validation study	9
3.1.4	Test results	9
3.1.5	Calculation and interpretation of relative trueness study	100
3.1.6	Conclusion (RT study)	15
3.2	Accuracy profile study	15
3.2.1	Categories, sample types and strains	16
3.2.2	Calculations and interpretation of accuracy profile study	17
3.3	Inclusivity / exclusivity	200
3.3.1	Protocols	200
3.3.2	Results	200
3.3.3	Conclusion	211
3.4	Limit of quantification (LOQ)	211
3.5	Conclusion (MCS)	211
4	Interlaboratory study	21
4.1	Study organisation	22
4.1.1	Collaborators	22
4.1.2	Matrix	22
4.1.3	Sample preparation	22
4.1.4	Labelling and shipping	22
4.1.5	Analysis of Samples	22
4.2	Experimental parameter controls	23
4.2.1	Detection of target organism in the matrix before inoculation	23
4.2.2	Sample variation	23
4.2.3	Logistic conditions	23
4.3	Calculation and summary of data	24
4.3.1	MicroVal Expert laboratory results	24
4.3.2	Results obtained by the collaborative laboratories	24
5	Overall conclusions of the validation study	27
	ANNEX A: Flow diagram of the reference and alternative method	27
	ANNEX B: Kit insert	27

1 Introduction

In this project a MicroVal validation study, based on ISO 16140-2:2016, of alternative method(s) for the enumeration of coagulase positive *Staphylococci* species in five different food categories was carried out by Campden BRI as the MicroVal Expert Laboratory.

This is a renewal of a method that was originally validated according to the superseded ISO16140:2003 standard for enumeration of coagulase-positive staphylococci in all foods. The original study was done by University of Ghent.

Five levels of contamination were used for the study, covering a minimum, a central and a maximum level plus two intermediary levels. Duplicate test portions were examined for each sample tested and this data can be partially used for the RT part of the renewal study as it covers two of the required food types within each category. Additional RT data was needed for a third type of food matrix within each category (Table 1).

There was no data available in the original study design to do the AP analysis as this part requires the testing of five replicate test portions and therefore new data was collected for this part (Table 4).

There is insufficient Incl/Excl data available from the original study to cover the requirements of ISO16140-2:2016 so some additional strains were tested in this renewal study. There were 36 inclusivity cultures tested so a further 14 were needed and there were 28 exclusivity cultures so a minimum of 2 more were required

The alternative method used was:

- Brilliance™ Staph 24 Agar for Enumeration of coagulase positive *Staphylococci* species following incubation at $37 \pm 1^\circ\text{C}$ for $24 \pm 2\text{h}$

The reference method used was:

- ISO 6888-1:1999 Microbiology of food and animal feeding stuffs- Horizontal method for the enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species) —Part 1: Technique using Baird-Parker agar
- ISO 6888-1:1999 DAM2:2017(E). Microbiology of food and animal feeding stuffs- Horizontal method for the enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species) —Part 1: Technique using Baird-Parker agar medium. AMENDMENT 2 Inclusion of an alternative confirmation procedure

Categories included:

- Dairy products
- Fishery products,
- Chocolate, bakery and confectionery products,
- Meat and meat products
- Multi-component foods

Criteria evaluated during the study have been:

- Relative trueness study;
- Accuracy profiles
- Inclusivity and exclusivity
- Interlaboratory Study

The final conclusion on the Method Comparison Study and ILS is summarized below:

The alternative method Brilliance™ Staph 24 Agar shows comparable performance to the reference methods (ISO 6888-1:1999, ISO 6888-1:1999 DAM2:2017(E) for the enumeration of coagulase positive *Staphylococci* species in a broad range of foods.

2 Method protocols

The Method Comparison Study was carried out using 10g gram portions of sample material.

According to ISO 16140-2 the reference method and alternative methods were performed with the same sample. The study was therefore a paired study design.

2.1 Reference method

See the flow diagram in Annex A.

Sample preparations used in the reference method were done according to ISO 6887-series parts 1, 2, 3, 4 and 5.

2.2 Alternative method

See the flow diagram of the alternative method in Annex A.

See the Brilliance™ Staph 24 Agar kit insert in Annex B.

The alternative method principle is based on chromogenic media:

Brilliance™ Staph 24 Agar: is a spread plate method intended to enumerate coagulase-positive Staphylococci. The agar is a transparent medium which is a highly selective and diagnostic chromogenic medium. Target organisms grow as dark blue colonies on a clear background. A picture is provided in Figure 1.

Figure 1: BSA



2.3 Study design

Samples of product containing the target organism were diluted 1 in 10 with an appropriate diluent according to ISO 6887 and homogenised in a stomacher.

Appropriate serial dilutions were made, and all relevant dilutions were analysed using the reference method and alternative method.

3 Method comparison study

3.1 Relative trueness study

The trueness study is a comparative study between the results obtained by the reference method and the results of the alternative method. This study was conducted using naturally or artificially contaminated samples. Different categories, types and items were tested for this.

A total of 5 categories were included in this validation study. A minimum of 15 items for each category were tested by both the reference method and the alternative method in the relative trueness study, with a minimum of 15 interpretable results per category.

Each category was made up of 3 types, with at least 5 items representative for each type.

3.1.1 Number of samples

The categories, the types and the number of samples analyzed are presented in Table 1.

Table 1 – Categories, types and number of samples analyzed

Category	Types		Number of interpretable samples			
			Old study		Renewal study	Total
Dairy products (combined category ; raw milk and heat processed)	a	Dry milk powders	10	0	4	14
	b	Raw milk and raw milk cheeses	10	6	0	16
	c	Dairy desserts e.g. chilled custard, trifle	0	0	5	5
	Total		20	6	9	35
RTE/RTRH Fishery products	a	Raw processed fish e.g. smoked mackerel	10	0	4	14
	b	RTE cooked fish i.e. prawns	10	0	4	14
	c	Frozen RTE / RTRH products	0	0	5	5
	Total		20	0	13	33
Chocolate, bakery products and confectionery	a	Confectionery: Cream patisserie, Chocolate mousse	10 10	0 0	4	24
	b	Dry powders e.g. cake mixes	0	0	5	5
	c	Low moisture products e.g. cakes, cookies, crackers	0	0	5	5
	Total		20	0	14	34
Meat (combined category; raw/RTC products and RTE/RTRH products)	a	Cured processed meat i.e. bacon	10	1	3	14
	b	Raw unprocessed meat i.e. mince	10	1	3	14
	c	RTE meats e.g. hams, pate	0	0	10	10
	Total		20	2	16	38
Multi component foods	a	RTE Deli salads with Mayonnaise	10	0	4	14
	b	RTRH chilled foods e.g. Pre-packed pancakes	10	0	4	14
	c	Composite processed meals	0	1	4	5
	Total		20	1	12	33
TOTAL DATA POINTS			100	9	64	173

3.1.2 Test sample preparation

It is preferable to test naturally contaminated samples. In order to attempt to use naturally contaminated samples, all fifteen samples from each category were first tested for the presence of naturally occurring target organism making a total of seventy five samples which were tested. In the original study there were 13 naturally positive samples for the RT study.

For the renewal study a further 71 samples were analysed for the RT. All these samples were screened for the presence of naturally occurring coagulase-positive staphylococci. However none of the sample were positive and so artificial contamination was needed. In order to ensure as wide a range of conditions were tested as possible, 15 different strains were used from 15 different food items and were either chilled frozen, heated or lyophilised.

All samples which were negative for natural contamination had <10cfu/g coagulase-positive staphylococci on both the reference method and alternative method.

Artificial contaminations were obtained by:

- Seeding with appropriate strains
 - o and storing chilled for minimum 48h at <5°C;
 - o and storing frozen for minimum 2 weeks at <-20°C or
 - o use of lyophilised cells, which were freeze dried, mixed into the dry powders and stored ambient for a minimum of 2 weeks before analysis

In the renewal study, the same strain was not used to inoculate more than 5 samples.

3.1.3 Protocols applied during the validation study

A single protocol was applied for the study.

Reference method plates were incubated at 37±1°C for 24±2h plus a further 24±2h.

Alternative method plates were incubated at 37±1°C for 24±2h.

In all cases the minimum incubation times were used.

Alternative method confirmation

Five typical colonies from each plate used in the calculations were confirmed using coagulase tests. Either ISO 6888-1:1999: BHI and rabbit plasma or ISO 6888-1:1999 DAM2:2017(E) - RPF Agar

3.1.4 Test results

The samples were analyzed by the reference and the alternative methods in order to have 15 interpretable results per incubation protocol, and 5 interpretable results per tested type by the two methods.



3.1.5 Calculation and interpretation of relative trueness study

The obtained data were analysed using the scatter plot. The graphs are provided with the line of identity ($y = x$).

Figures 2 to 7 shows the scatter plots for the individual categories and all categories.

Figure 2 - Scatter plot of the reference method versus alternative method results for Milk and dairy products




Figure 3- Scatter plot of the reference method versus alternative method results for Chocolate and bakery




Figure 4- Scatter plot of the reference method versus alternative method results for Meat (combined category; raw/RTC products and RTE/RTRH products)




Figure 5- Scatter plot of the reference method versus alternative method results for Ready to eat fish




Figure 6- Scatter plot of the reference method versus alternative method results for Multi component foods





Figure 7 - Scatter plot of the reference method versus alternative method results for all categories



According to ISO 16140-2:2016 6.1.2.3 the results of the scatter plot are interpreted based on a visual observation on the amount of bias and extreme results.

There was extremely good agreement between the two methods with almost no positive or negative bias.

A summary of the calculated values per category is provided in Table 2 and the Bland-Altman difference plot

Table 2- Summary of the calculated values per category –combined renewal and original data

Category.	n	\bar{D}	s_D	95% Lower limit	95% Upper limit
chocolate and bakery	34	-0.038	0.119	-0.283	0.208
Dairy	35	-0.025	0.128	-0.289	0.239
Meat	38	0.011	0.085	-0.164	0.187
multi component foods	35	0.039	0.134	-0.239	0.317
RTE/RTRH fish	33	-0.012	0.127	-0.274	0.250
All Categories	173	-0.005	0.121	-0.244	0.234

\bar{D} : Average difference SD: standard deviation of differences n: number of samples

Figure 8 – Bland-Altman difference plot for all the samples




Table 3 - Data which are outside of the accepted limits

Category	Types	Code	Food item	strain	Spiking/ seeding	Log (Ref)	Log (Alt)	Mean	Difference
RTE/RTRH fish	raw processed fish	10	smoked mackerel	<i>S.aureus</i> 1208	chill 2-3 days	2.114	2.653	2.384	0.539
chocolate and bakery	Dry powders	30	sponge mix	<i>S.aureus</i> ATCC 29213	Lyophilised cells	3.556	3.230	3.393	-0.326
chocolate and bakery	Confection ery	23	cream cake	<i>S.aureus</i> 1215	chill 2-3 days	5.898	5.544	5.721	-0.354
chocolate and bakery	Dry powders	29	muffin mix	<i>S.aureus</i> ATCC 29213	Lyophilised cells	3.322	3.000	3.161	-0.322
multi component foods	RTE deli salads	48	chargrilled sweetcorn coleslaw	<i>S.aureus</i> 3098	chill 2-3 days	3.886	4.146	4.016	0.260
multi component foods	RTE deli salads	47	salad with sweetcorn	<i>S.aureus</i> 3098	chill 2-3 days	2.146	2.695	2.420	0.548
multi component foods	RTRH chilled foods	54	southern fried chicken wings	<i>S.aureus</i> 1994	chill 2-3 days	5.613	5.301	5.457	-0.312
Dairy	raw cheese	170	Raw milk	Natural strain	chill 2-3 days	1.146	1.519	1.332	0.372
Dairy	raw cheese	174	Raw milk cheese	Natural strain	chill 2-3 days	4.204	3.708	3.956	-0.497

It is expected that not more than one in 20 data values will lie outside the CLs.

In this study there were 9 data points from a total of 173 data points which were outside of the accepted limits. This meets the expectation. The data covered 4 different food categories, 5 different *Staphylococcus aureus* strains, and naturally present coagulase positive staphylococci. The differences between the results of the compared methods are in most of the case below the absolute value of 0,5 Log CFU/g, and only two are slightly above.

3.1.6 Conclusion (RT study)

The relative trueness of the Alternative method is satisfied as the expectation of not more than 1 in 20 data points outside of the acceptability limits is met

3.2 Accuracy profile study

The accuracy profile study is a comparative study between the results obtained by the reference and the results of the alternative method. This study is conducted using artificially contaminated samples, using one type per category.

3.2.1 Categories, sample types and strains

It is possible to run this study in two different ways. It is possible to use either 2 separate batches of a single item for each food type. Or it is possible to use a single batch of 2 different items for each food type.

In this study five food categories were tested with a single batch of two different food types using 6 samples per type. Two samples were contaminated at a low level, 2 at intermediate level, 2 at a high level. For each sample, 5 replicates (5 different test portions) were tested. A total of 30 samples were analysed per food type. The conditions tested are shown in Table 4

Each sample was bulk inoculated and five replicate test portions examined from the bulk sample.

Table 4 - Categories, types, items, strains and inoculation levels for accuracy profile study

Category	Types	Strain	Item	Level
Dairy products (combined category; raw milk and heat processed)	Dairy desserts	<i>S. aureus</i> CRA 1215 from cheese	Chilled custard	Level 1x5: 100-250 cfu /g
				Level 2x5: 1000- 5000 cfu/g
				Level 3x5 10,0000- 50,0000 cfu/g
			Whipped cream	Level 1x5: 100-250 cfu /g
				Level 2x5: 1000- 5000 cfu/g
				Level 3x5 10,0000- 50,0000 cfu/g
RTE/RTRH Fishery products	RTE fish products	<i>S. aureus</i> CRA 1208 from smoked fish	Smoked salmon	Level 1x5: 100-250 cfu /g
				Level 2x5: 1000- 5000 cfu/g
				Level 3x5 10,0000- 50,0000 cfu/g
			Tuna pate	Level 1x5: 100-250 cfu /g
				Level 2x5: 1000- 5000 cfu/g
				Level 3x5 10,0000- 50,0000 cfu/g
Chocolate, bakery products and confectionery	Pastries	<i>S. aureus</i> CRA 2078 Milk powder	Chilled patisserie	Level 1x5: 100-250 cfu /g
				Level 2x5: 1000- 5000 cfu/g
				Level 3x5 10,0000- 50,0000 cfu/g
			Chocolate filled pain au chocolate	Level 1x5: 100-250 cfu /g
				Level 2x5: 1000- 5000 cfu/g
				Level 3x5 10,0000- 50,0000 cfu/g
Meat (combined category; raw/RTC products and RTE/RTRH products)	RTE meats	<i>S. aureus</i> CRA 1219 beef	Sliced ham	Level 1x5: 100-250 cfu /g
				Level 2x5: 1000- 5000 cfu/g
				Level 3x5 10,0000- 50,0000 cfu/g
			Pastrami	Level 1x5: 100-250 cfu /g
				Level 2x5: 1000- 5000 cfu/g
				Level 3x5 10,0000- 50,0000 cfu/g
Multi component foods	Composite foods with raw /processed ingredients	<i>S. aureus</i> CRA 5932 from pasta	Pasta salad	Level 1x5: 100-250 cfu /g
				Level 2x5: 1000- 5000 cfu/g
				Level 3x5 10,0000- 50,0000 cfu/g
			Sandwich spread	Level 1x5: 100-250 cfu /g
				Level 2x5: 1000- 5000 cfu/g
				Level 3x5 10,0000- 50,0000 cfu/g

*these are target values only and actual values may be ± 1 log from the target dependent on microbial behaviour

Total number of samples tested= 150

3.2.2 Calculations and interpretation of accuracy profile study

The statistical results and the accuracy profiles are provided in Figures 9 to 13.

The calculations were done using the AP Calculation Tool MCS (Clause 6-1-3-3 calculation and interpretation of accuracy profile study) available on <http://standards.iso.org/iso/16140>

Figure 9 Accuracy profile for Category : Milk and dairy products (types custard and cream)




Figure 9 Accuracy profile for Category : Chocolate and bakery(types patisserie)




Figure 10 Accuracy profile for Category : Raw poultry and meats (types RTE meats)




Figure 11 Accuracy profile for Category : Ready to eat fish (types RTE fish)





Figure 12 Accuracy profile for Category : Multi component foods (types products with mayonaïse)



If any of the upper or lower limits exceeded the 0.5log AP limits and the standard deviation of the reference method was >0.125, additional evaluation procedures are required, as described in ISO 16140-2:2016 and the new acceptability limits are calculated

In this study all four categories met the AL of 0.5log. These were; Dairy products; Fishery products; Chocolate, bakery and confectionery products; Meat and meat products.

For Multi-component foods, the standard deviation of the reference method was >0.125 and one out of 12 limits was exceeded. This was the upper CL for the low level pasta salad.

The additional calculations were carried out and the reference method met the newly calculated AL of ± 0.684 .

The accuracy of the Alternative method is satisfied as the all categories met the 0.5log AL or in the case of Multi-component foods, the re-calculated AL.

3.3 Inclusivity / exclusivity

Inclusivity is the ability of the alternative method to detect the target analyte from a wide range of strains.

Exclusivity is the lack of interference from a relevant range of non-target strains of the alternative method.

In this study, the target organism is coagulase positive Staphylococci which consists of three species:- *Staphylococcus aureus*, - *Staphylococcus hyicus*, - *Staphylococcus intermedius*.

3.3.1 Protocols

- **Inclusivity**

A minimum of fifty strains of coagulase positive staphylococci needed to be tested to meet the requirements of ISO16140-2:2016.

Thirty four strains of coagulase positive staphylococci were evaluated in the original study.

In addition, a further seventeen cultures were tested here making a total of fifty one strains tested.

In this renewal study each strain was grown in Nutrient Broth in at $30\pm 1^\circ\text{C}$ for 18-24h and appropriate dilutions were made for testing. Each strain was tested once with the alternative method, the reference method and a non-selective agar.

- **Exclusivity**

A minimum of thirty strains of non target organisms needed to be tested to meet the requirements of ISO16140-2:2016. Twenty eight strains of non-coagulase positive staphylococci were evaluated in the original study.

In addition, a further five cultures were tested here making a total of thirty three strains tested. All strains were grown in appropriate selective broths and incubation conditions and appropriate dilutions were made for testing. Each strain was tested once with the alternative method, the reference method and a non-selective agar.

3.3.2 Results

- Inclusivity

Of the 51 inclusivity strains tested, 48 were detected by the reference method and the alternative method. One strain (*Staph. intermedius* code 407) was detected by the alternative method but not the reference method. And two strains (*Staph. hyicus* code 285 , *Staph. hyicus* code 406) were not detected by either method.

- Exclusivity

Of the 35 exclusivity strains tested, 2 were detected by the reference method and the alternative method. (*Staph. caprae* code 21 from Goats Milk; *Staph. epidermidis* code 402 from Goats Milk).

A further three strains were detected by the reference method only (*Staph. carnosus*, CRA 284 from Fermented sausage; *Staph saprophyticus*, CRA 3191 from dry sausage; *Staph. warneri* CRA 262 from Salami) but these colonies were not confirmed and therefore would not be detected as coagulase positive Staphylococci

3.3.3 Conclusion

The alternative method Brilliance™ Staph 24 Agar for enumeration of coagulase positive *Staphylococci* species in foods was shown to be specific and selective and give comparable performance to the reference method.

3.4 Limit of quantification (LOQ)

The limit of Quantification (LOQ) is only required for instrumental measurements. It was not done in this study

3.5 Conclusion (MCS)

Overall, the conclusions for the Method Comparison are:

- The alternative method Brilliance™ Staph 24 Agar for enumeration of coagulase positive *Staphylococci* species shows satisfactory results for relative trueness;
- The alternative method Brilliance™ Staph 24 Agar for enumeration of coagulase positive *Staphylococci* species shows satisfactory results for accuracy profile;
- The alternative method Brilliance™ Staph 24 Agar for enumeration of coagulase positive *Staphylococci* species is selective and specific.

4 Interlaboratory study

The inter-laboratory study is a study performed by multiple laboratories testing identical samples at the same time, the results of which are used to estimate alternative-method performance parameters.

4.1 Study organisation

4.1.1 Collaborators

Samples were sent to 11 laboratories in four different countries.

4.1.2 Matrix and strain used

Pasteurised milk samples were inoculated with *Staph. aureus* LFMFP-UGent N°532 isolated from raw milk cheese. Samples were individually inoculated with the relevant dilution of the *Staph. aureus* strain.

4.1.3 Sample preparation

Samples were prepared and inoculated and despatched as described below:

Each laboratory received eight samples of 25 g, i.e. two samples per inoculation level L1, L2, and L3 and 2 samples of the blank L0. All samples were blind-coded.

The target contamination levels were:

Level 0 < 10 cfu/g (Blank)

Level 1: 500 cfu/ml

Level 2: 5000 cfu/ml

Level 3: 50,000cfu/ml

A set of samples was also prepared for the EL although the data from these was not used in the data analysis.

4.1.4 Labelling and shipping

Before dispatch, each set of samples were packed into isothermal boxes, containing cooling blocks, and express-shipped to the different laboratories. It was aimed to have the samples delivered within 24-48h.

Also, an additional vial containing a temperature monitoring sensor was packed with each set of samples. These were sent back to the expert lab to provide a profile of the temperature during transportation and receipt at the collaborators.

Upon receipt, each collaborating laboratory tested each sample according to the reference method and alternative method following the instruction sheets which were sent prior to the start of the study. In addition, the organising laboratory tested a set of samples at the same time as the collaborating laboratories to confirm the presence of the target organisms and the contamination levels. The expert laboratory data was not used in the calculations.

4.1.5 Analysis of Samples

Collaborative study laboratories and the expert laboratory carried out the analyses on Wednesday October 7th, 2009 with the alternative and reference methods. The analyses by the reference method and the alternative method were performed on the same day.

4.2 Experimental parameter controls

4.2.1 Detection of coagulase positive staphylococci in the matrix before inoculation

In order to ensure the absence of coagulase positive staphylococci in the food matrix, the reference method was performed on five portions (25 g) before the inoculation. All the results were negative.

4.2.2 Strain stability during transport

Duplicate samples of each of L1, L2 and L3 were tested at the EL for stability after storage at <8°C for 0, 24 and 48h on both the alternative and reference method.

Table 5: Levels of coagulase positive staphylococci (\log_{10} cfu/g) in stability samples stored at 4°C

Level and time	Reference: BPA			Alternative: BSA		
	0h	24h	48h	0h	24h	48h
low a	2.51	2.68	2.56	2.56	2.68	2.64
low b	2.54	2.66	2.61	2.56	2.72	2.66
medium a	3.62	3.48	3.53	3.73	3.47	3.61
medium b	3.58	3.53	3.51	3.53	3.51	3.69
high a	4.48	4.50	4.51	4.48	4.50	4.53
high b	4.45	4.43	4.54	4.53	4.42	4.36

The data showed that the levels of coagulase positive staphylococci were not affected by the storage conditions and were stable during chill storage with no increase after 48h at 4°C.

4.2.3 Logistic conditions

The temperatures measured at receipt by the collaborators, the temperatures registered by the thermo-probe, and the receipt dates are given in Table 6.

Table 6 - Sample temperatures at receipt

Organising laboratory	Temperature measured at receipt (°C)	Receipt date and time	Analysis date
1	5.0°C	Day 0 11.10	7 th October 2009
2	3.5°C	Day 0 10.30	7 th October 2009
3	7.3°C	Day 0 12.00	7 th October 2009
4	3.8°C	Day 1 12.15	7 th October 2009
5	4.9°C	Day 0 09.45	7 th October 2009
6	6.6°C	Day 1 11.00	7 th October 2009
7	9.5°C	Day 0 13.00	7 th October 2009
8	Forgotten to measure	Day 0 11.15	7 th October 2009
9	5.0°C	Day 1 11.15	7 th October 2009
10	1.6°C	Day 0 11.00	7 th October 2009
11	3.4°C	Day 1 11.00	7 th October 2009
EL	5.0°C	Day 0 08.00	7 th October 2009

No problems were encountered during the transport or at receipt for the 11 collaborators.

All the samples were delivered on time and in appropriate conditions.

Temperatures during shipment and at receipt were all correct. The temperature reading at receipt from the water sample was <8°C except for laboratory 7 which had a measured temperature of 9.5°C

4.3 Calculation and summary of data

4.3.1 MicroVal Expert laboratory results

The results obtained by the expert laboratory are given in Table 7.

Table 7 – Results obtained by the expert lab(cfu/g)

Level	Reference method	Alternative method
Blank	<10	<10
Low	2.67	2.68
Low	2.66	2.72
Medium	3.48	3.48
Medium	3.53	3.52
High	4.51	4.51
High	4.43	4.41

4.3.2 Results obtained by the collaborative laboratories

The data from the collaborative trial were calculated and interpreted according to section 6.2.3 of ISO 16140-2:2016 using the freely available Excel® spreadsheet (<http://standards.iso.org/iso/16140>). Version 14-03-2016 was used for these calculations.

The results obtained by the collaborators are shown in Table 8.


The accuracy profile plot is shown in Figure 13 and the statistical analysis of the data shown in Table 9.

Table 8: Summary of the results of the interlaboratory study per analyte level

Collaborator/level		Reference method (Log cfu/g)		Alternative method (Log cfu/g)	
		Duplicate 1	Duplicate 2	Duplicate 1	Duplicate 2
1	low	2.60	2.57	2.53	2.62
2	low	2.80	2.81	2.84	2.83
3	low	2.85	2.85	2.85	2.83
4	low	2.72	2.74	2.76	2.72
5	low	2.52	2.58	2.60	2.52
6	low	2.92	2.82	2.76	2.86
7	low	2.65	2.41	2.65	2.48
8	low	2.40	2.43	2.51	2.49
9	low	2.45	2.48	2.45	2.48
10	low	2.54	2.51	2.59	2.57
11	low	2.60	2.72	2.70	2.59
1	medium	3.60	3.81	3.72	3.86
2	medium	3.91	3.72	3.63	3.73
3	medium	3.80	3.85	3.85	3.90

4	medium	3.80	3.89	3.83	3.81
5	medium	3.52	3.51	3.54	3.40
6	medium	3.79	3.84	3.81	3.81
7	medium	3.59	3.69	3.64	3.53
8	medium	3.45	3.41	3.45	3.51
9	medium	3.32	3.49	3.48	3.45
10	medium	3.62	3.45	3.51	3.49
11	medium	3.49	3.45	3.56	3.52
1	high	4.56	4.86	4.68	4.88
2	high	4.95	4.65	4.91	4.72
3	high	4.88	4.70	4.81	4.75
4	high	4.90	4.85	4.91	4.88
5	high	4.56	4.54	4.58	4.45
6	high	4.82	4.95	4.85	4.93
7	high	4.66	4.69	4.69	4.54
8	high	4.56	4.46	4.58	4.40
9	high	4.36	4.48	4.51	4.41
10	high	4.71	4.72	4.72	4.69
11	high	4.53	4.59	4.48	4.46

Figure 14. Accuracy profile of Brilliance™ Staph 24 Agar from the ILS



**Standardized report - Quantitative methods -
Method Comparison Study and ILS
2008LR11 renewal Brilliance™ Staph 24 Agar
Summary report**



Table 9. Statistical analysis of the ILS data according to the ISO spreadsheet

Accuracy profile

Study Name

Date

Coordinator

Tolerance probability (beta)

Acceptability limit in log (lambda)

Staph Brilliance

calculated 09/01/2018 from existing data

Campden BRI

80%80%80%

0.500.500.50

Alternative method

Levels

Low

Medium

High

Target value

Number of participants (K)

Average for alternative method

Repeatability standard deviation (sr)

Between-labs standard deviation (sL)

Reproducibility standard deviation (sR)

Corrected number of dof

Coverage factor

Interpolated Student t

Tolerance interval standard deviation

Lower TI limit

Upper TI limit

Bias

Relative Lower TI limit (beta = 80%)

Relative Upper TI limit (beta = 80%)

Lower Acceptability Limit

Upper Acceptability Limit

2.635

3.636

4.680

11

2.646

3.637

4.674

0.057

0.057

0.086

0.130

0.156

0.158

0.142

0.166

0.180

11.746

11.237

12.570

1.414

1.419

1.406

1.358

1.362

1.353

0.1477

0.1734

0.1874

2.445

3.401

4.420

2.847

3.873

4.927

0.011

0.001

-0.007

-0.190

-0.235

-0.260

0.212

0.237

0.247

-0.50

-0.50

-0.50

0.50

0.50

0.50

New acceptability limits may be based on reference method pooled variance

Pooled repro standard dev of reference

0.173

FALSE

Application of clause 6.2.3

Step 8: If any of the values for the β -ETI fall outside the acceptability limits, calculate the pooled average reproducibility standard deviation of the reference method.

Step 9: Calculate new acceptability limits as a function of this standard deviation.

Reference method

Low

Medium

High

11

11

11

2.635

3.636

4.680

0.063

0.086

0.109

0.152

0.161

0.132

0.164

0.182

0.172

11.604

12.489

14.845

FALSE

FALSE

Select ALL blue lines to draw the accuracy profile as illustrated in the worksheet "Graph Profile"

5 Overall conclusions of the validation study

- The alternative method Brilliance™ Staph 24 Agar for enumeration of coagulase positive Staphylococci species shows satisfactory results for relative trueness;
- The alternative Brilliance™ Staph 24 Agar for enumeration of coagulase positive Staphylococci species shows satisfactory results for accuracy profile;
- The alternative Brilliance™ Staph 24 Agar for enumeration of coagulase positive Staphylococci species is selective and specific.
- The alternative Brilliance™ Staph 24 Agar for enumeration of coagulase positive Staphylococci species shows satisfactory performance in the ILS

The alternative Brilliance™ Staph 24 Agar shows comparable performance to the reference method ISO 6888-1:1999 DAM2:2017 for enumeration of coagulase positive Staphylococci species

Date, 28/03/2019


Signature



Annexes

- A. Flow diagram of the reference and alternative method
- B. Test kit insert

ANNEX A: Flow diagram of the alternative method and reference methods



Calculate cfu/g taking into account the number of confirmed positive colonies

*It is possible to spread 1ml of the initial suspension on 3 plates (90mm) for low number estimation. This should be done in duplicate

Note RPF agar and Tube coagulase both used for inclusivity exclusivity cultures. RPF only used for RT and AP

**Standardized report - Quantitative methods -
Method Comparison Study and ILS
2008LR11 renewal Brilliance™ Staph 24 Agar
Summary report**



ANNEX B: Kit insert(s) -latest version provided as a separate document

Oxoid Prepared Medium Brilliance™ Staph 24 Agar

REF: PO1186A

Intended Use
IVD

Brilliance Staph 24 is a selective, chromogenic medium for the isolation and enumeration of coagulase positive staphylococci (CPS) in clinical and food samples. For professional use only (in vitro diagnostic use).

Summary and Explanation
Brilliance Staph 24 Agar is a selective and diagnostic chromogenic medium for the isolation and enumeration of coagulase-positive staphylococci (CPS) in foods and clinical samples, within 24 hours.
With traditional media for food, such as Baird-Parker Egg Yolk Tellurite Agar, it takes up to 48 hours to obtain a result.
With traditional media for clinical samples, such as Columbia CNA Agar, detection is not limited to staphylococci, so organisms like streptococci can also grow.
Brilliance Staph 24 reduces non-target organism growth while allowing all strains of CPS to grow uninhibited, providing more accurate results and reducing the number of confirmatory tests required.

Principle
Brilliance Staph 24 Agar provides dark blue colonies on a clear agar background after only 24 hours incubation. This contrast makes identification and enumeration of CPS simple. Brilliance Staph 24 Agar is designed to detect all species of CPS, not just *Staphylococcus aureus*, to make this agar an all-inclusive test. The design also restricts non-target organism growth while allowing all strains of CPS to grow uninhibited, leading to more accurate enumeration and a reduction in the number of confirmation tests required.
Sodium pyruvate is added to improve the recovery of stressed cells, and a carefully formulated blend of peptones and growth factors ensures rapid growth of target organisms. The chromogen is specifically activated by CPS which colours positive colonies dark blue, while coagulase-negative staphylococci are inhibited or remain colourless.

Typical Formula*

	grams per litre
Peptone	25.0
Chromogenic mix	5.0
Sodium pyruvate	4.0
Lithium chloride	5.0
Antibiotic mix	20.0 ml
Agar	14.0

* Adjusted as required to meet performance standards

Physical Characteristics
Colour: Straw
Clarity: Clear
Fill weight: 17.0
pH: 7.2 ± 0.2

Precautions
This product is for *in vitro* diagnostic use and should only be used by trained individuals. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products. It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully.
For professional use only. Safety Data Sheet available on request.

Storage
This product is ready to use and no further preparation is necessary.
Store product in its original packaging at 6–12°C until used. Store away from light.
Allow product to equilibrate to room temperature before use. Do not incubate prior to use.

Specimen Collection, Handling and Storage
Specimen should be collected and handled following the recommended guidelines.

Materials Required but Not Supplied
(1) Inoculating loops, swabs, collection containers
(2) Incubators
(3) Quality control organisms
More information on www.thermoscientific.com/microbiology

Procedure
Food application
(1) Dilute sample in an appropriate diluent.
Spread 0.1 mL aliquots over the surface of the Brilliance Staph 24 plate until absorbed.....
(2) Incubate for 24h ± 2h at 37°C ± 1°C.
(4) Blue colonies are presumptive positive for coagulase positive staphylococci
(5) Identifications as *S. aureus* can be confirmed coagulase tests

Clinical application
Brilliance Staph 24 Agar can be inoculated from screening swabs taken from hospital patients or staff, isolated colonies or from liquid suspension such as *Contras™* MRSA broth (EB1225B), according to local guidelines.
(1) Inoculate the Brilliance Staph 24 Agar etc.
(2) Incubate plates aerobically for 18–24 hours at 36 ± 1°C.
(4) Blue colonies are presumptive positive for coagulase positive staphylococci
(5) Identifications as *S. aureus* can be confirmed with Staphytest Plus or Dryspot Staphytest Plus
More information on: www.thermofisher.com

Quality Control
Results after incubation at 35–39°C for 18–24 hours

Control strain	Growth characteristic
----------------	-----------------------

**Standardized report - Quantitative methods -
Method Comparison Study and ILS
2008LR11 renewal Brilliance™ Staph 24 Agar
Summary report**



<i>Staphylococcus aureus</i> ATCC® 25923	Blue colonies
<i>Staphylococcus aureus</i> ATCC® 6538	Blue colonies
<i>Staphylococcus saprophyticus</i> ATCC® 15305	White colonies
<i>Escherichia coli</i> ATCC® 25922	No growth
<i>Bacillus cereus</i> ATCC® 10876	No growth
<i>Enterococcus faecalis</i> ATCC® 29212	No growth
<i>Staphylococcus epidermidis</i> ATCC® 12228	Inhibited as slight background film

Note:

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The product should not be used if:

- (1) The product is contaminated
- (2) The colour has changed
- (3) The expiration date has passed
- (4) There are other signs of deterioration

Performance

MicroVal ISO 16140 study*

Linearity study:

A total of 120 samples (consisting of meat products, fish products, bakery products, composite food and dairy products) were artificially contaminated. Each food type was contaminated within the ranges of: 1×10^2 to 5×10^2 CFU/g, 5×10^2 to 5×10^3 CFU/g, 5×10^3 to 5×10^4 CFU/g, 5×10^4 to 5×10^5 CFU/g and 5×10^5 to 5×10^6 CFU/g. Samples were analysed by both the ISO 6888-1:1999/Amd 1:2003 method and Brilliance Staph 24 Agar method.

Food category	Food product/strain	Regression line
Meat	Bacon	$y = 0.086 + 1.033 x$
	Mixed meat	$y = 0.013 + 1.006 x$
Fish products	Cooked peeled shrimp	$y = 0.080 + 0.977 x$
	Smoked mackerel	$y = 0.120 + 0.978 x$
Bakery products	Chocolate mousse	$y = 0.077 + 0.983 x$
	Pasty cream	$y = 0.045 + 0.956 x$
Composite food	Pre-packed pancakes	$y = 0.104 + 0.975 x$
	Raw mayonnaise based salad	$y = 0.061 + 0.996 x$
Dairy	Dried milk powder	$y = 0.030 + 0.987 x$
	Unpasteurised milk based cheese	$y = 0.045 + 0.968 x$

Relative accuracy study

One hundred and twenty artificially contaminated samples and 18 naturally contaminated samples were tested. The contamination (concentration) ranges were as follows:

Food category	Contamination range (in log CFU/g)
Meat products	2.30-6.57

Fish products	2.16-6.32
Bakery products	2.60-6.29
Composite food	2.30-6.54
Dairy products	2.64-6.58

The equation of the regression line between the Brilliance Staph 24 Agar and the ISO reference method for all categories combined is:

$$y = 0.0074 + 0.9981 x$$

$$R^2 = 0.9985$$

$$y = \log(N \text{ Brilliance Staph 24 Agar})$$

$$x = \log(N \text{ ISO reference method})$$

The linearity and relative accuracy study proved Brilliance Staph 24 Agar to have a good degree of equivalence to Baird Parker Egg Yolk Tellurite Agar medium for the enumeration of coagulase positive staphylococci.

Selectivity (inclusivity/exclusivity)

A total of 38 coagulase positive staphylococci (31 *S. aureus*, 2 *S. intermedius*, 2 *S. caprae*, 1 *S. epidermidis*) were tested, all gave positive results on Brilliance Staph 24 Agar. Twenty eight strains of coagulase negative staphylococci and non-staphylococci gave a negative result with Brilliance Staph 24 Agar. In comparison, 13 of the exclusivity isolates produced typical colonies on Baird Parker Egg Yolk Tellurite Agar medium which required confirmation. These were later identified by the tube-coagulase test as coagulase negative.

Practicability

The Brilliance Staph 24 Agar medium is an agar based culture medium shown to be of practical use in the lab as the blue typical colonies enabled convenient reading of the plates and enumeration of presumptive coagulase positive staphylococci.

Interlaboratory study

Eleven different laboratories tested samples of pasteurised milk that were artificially contaminated with a single strain to provide samples with low (5×10^2 CFU/ml), medium (5×10^3 CFU/ml) and high (5×10^5 CFU/ml) contamination levels. Uninoculated samples were used to provide the fourth contamination level (0 CFU/ml). Each laboratory received duplicate blind-coded samples for each contamination level which were tested by both the ISO 6888-1:1999/Amd 1:2003 method and Brilliance Staph 24 Agar method.

Contamination level	No. samples taken into account	ISO reference method		Brilliance Staph 24 Agar		Bias
		r (sd)	R (sd)	r (sd)	R (sd)	
Low	22	0.031	0.189	0.042	0.146	0.019
Medium	22	0.094	0.279	0.052	0.197	-0.005
High	22	0.105	0.202	0.105	0.281	-0.010

r (sd) = repeatability
R (sd) = reproducibility

The results from the method comparison study and the interlaboratory study revealed that there were no substantial differences between Brilliance Staph 24 Agar and Baird Parker Egg Yolk Tellurite Agar medium according to the reference method ISO 6888-1:1999/Amd 1:2003 for enumeration of coagulase positive Staphylococci.

Clinical study

Five hundred fifty-eight wound swabs from a UK hospital were tested on Brilliance Staph 24 Agar and a competitor chromogenic agar for detection of *S. aureus* (plate 'A'). Performance of the two media is summarized below:

Performance	Brilliance	Plate 'A'
-------------	------------	-----------

Thermo
SCIENTIFIC
A Thermo Fisher Scientific Brand