

2022LR113 Neogen Soleris Non-Fermenting Total Viable Count (NF-TVC) vial MCS ILS Summary Report 14/10/2024



**Method Comparison Study Report for the ISO 16140-2:2016 validation of Neogen Soleris Non-Fermenting Total Viable Count (NF-TVC), for the detection of Total Viable Count at a threshold of 10 cfu per g in a broad range of foods**

MicroVal study number: 2022LR113

Method/Kit name: Neogen Soleris® Non-Fermenting Total Viable Count (NF-TVC)

Report version: MCS ILS Summary Report 14/10/2024

MicroVal Expert Laboratory:

Suzanne Jordan and Alice Foxall  
Station Road,  
Chipping Campden,  
Gloucs,  
GL55 6LD, UK  
Tel: 0044 1386 842000  
Email: [suzanne.jordan@campdenbri.co.uk](mailto:suzanne.jordan@campdenbri.co.uk)  
[www.campdenbri.co.uk](http://www.campdenbri.co.uk)

## Foreword

This report is prepared in accordance with ISO 16140-2:2016 and MicroVal technical committee interpretation of ISO 16140-2 v2.5

Company: Neogen Corporation  
  
The Dairy School,  
Auchincruive,  
Ayr,  
KA6 5HU, Scotland, UK

Expert Laboratory: Campden BRI  
Station Road  
Chipping Campden  
Gloucs,  
GL55 6LD, UK

Method/Kit name: Neogen Soleris® Non-Fermenting Total Viable Count (NF-TVC)

Validation standard: Microbiology of the food chain— Method validation

- Part 1: Vocabulary (ISO 16140-1:2016)
- Part 2: Protocol for the validation of alternative (proprietary) methods against a reference method (ISO 16140-2:2016)

Reference method: ISO 4833-1:2013 Microbiology of the food chain - Horizontal method for the enumeration of microorganisms - Part 1: Colony count at 30°C by the pour plate technique

Scope of validation: Broad range of foods

Certification organization: LRQA

### List of abbreviations

|            |   |
|------------|---|
| A(It)      | Alternative method  |
| AL         | Acceptability Limit   |
| Art. Cont. | artificial contamination  |
| CFU        | Colony Forming Units  |
| EL         | Expert Laboratory   |
| FP         | False Positive  |
| FPR        | False Positive Ratio  |
| g          | Gram  |
| h          | Hour  |
| ILS        | Interlaboratory Study   |
| LOD        | Level of Detection  |
| MCS        | Method Comparison Study   |
| min        | minute  |
| ml         | millilitre  |
| MR         | (MicroVal) Method Reviewer  |
| MVTC       | MicroVal Technical Committee  |
| NA         | Negative Agreement  |
| na         | not applicable  |
| ND         | Negative Deviation  |
| neg (-)    | negative/no growth/no reaction/target not detected                              |
| NS         | Non-Suspect growth  |
| nt         | not tested  |
| PA         | Positive Agreement  |
| PD         | Positive Deviation  |
| pos (+)    | positive/growth/target detected   |
| PPNA       | Presumptive Positive Negative Agreement (belongs to the False Positive results) |
| PPND       | Presumptive Positive Negative Deviation (belongs to the False Positive results) |
| R(ef)      | Reference method  |
| RLOD       | Relative Level of Detection   |
| RT         | Relative Trueness   |
| S          | Suspect growth  |
| SE         | Relative Sensitivity  |
| SP         | Relative Specificity  |
| TP         | True Positive   |

## Contents

|          |  |           |
|----------|--|-----------|
| <b>1</b> | <b>Introduction</b>                                  | <b>7</b>  |
| <b>2</b> | <b>Method protocols</b>                              | <b>9</b>  |
| 2.1      | Reference method                                     | 9         |
| 2.2      | Alternative method                                   | 9         |
| 2.3      | Study design   | 9         |
| <b>3</b> | <b>Method comparison study</b>                       | <b>10</b> |
| 3.1      | Sensitivity Study                                    | 10        |
| 3.1.1    | Categories and sample types                          | 10        |
| 3.1.2    | Test sample preparation                              | 12        |
| 3.1.3    | Confirmation protocols                               | 12        |
| 3.1.4    | Sensitivity study results                            | 13        |
| 3.1.5    | Sensitivity study calculations                       | 14        |
| 3.1.6    | Discordant results                                   | 16        |
| 3.1.7    | Conclusion sensitivity study                         | 19        |
| 3.2      | Relative level of detection study                    | 19        |
| 3.2.1    | Categories, sample types and strains                 | 19        |
| 3.2.2    | Test sample preparations                             | 20        |
| 3.2.3    | RLOD study results                                   | 20        |
|          | The LOD50 per Category is given in Tables 12 and 13. | 20        |
| 3.2.4    | Conclusion RLOD study                                | 21        |
| 3.3      | Inclusivity/exclusivity study                        | 21        |
| 3.3.1    | Protocols  | 21        |
| 3.3.2    | Results inclusivity and exclusivity study            | 22        |
| 3.3.3    | Conclusion inclusivity and exclusivity study         | 22        |

|          |  |           |
|----------|--|-----------|
| <b>4</b> | <b>Conclusions Method Comparison Study</b>   | <b>22</b> |
| <b>5</b> | <b>Interlaboratory study</b>   | <b>23</b> |
| 5.1      | Study organisation   | 23        |
| 5.2      | Experimental parameters controls   | 24        |
| 5.2.1    | Detection of total viable count in the matrix before inoculation   | 24        |
| 5.2.2    | Strain stability during transport  | 24        |
| 5.2.3    | Logistic conditions  | 25        |
| 5.3      | Calculation and summary of data  | 27        |
| 5.3.1    | MicroVal Expert laboratory results   | 27        |
| 5.3.2    | Results obtained by the collaborative laboratories   | 27        |
| 5.3.3    | Results of the collaborators retained for interpretation   | 30        |
| 5.3.4    | Calculation of the specificity percentage (SP)   | 32        |
| 5.3.5    | Calculation of the sensitivity ( $SE_{alt}$ ), the sensitivity for the reference method ( $SE_{ref}$ ), the relative trueness (RT) and the false positive ratio for the alternative method (FPR) | 32        |
| 5.3.6    | Interpretation of data   | 33        |
| 5.3.7    | Evaluation of the RLOD between laboratories  | 35        |
| 5.4      | Root cause analysis – interlaboratory study  | 35        |
| 5.4.1    | Recruitment of Participants  | 36        |
| 5.4.2    | Positive blank samples   | 38        |
| 5.4.3    | Summary  | 44        |
| <b>6</b> | <b>CONCLUSION</b>  | <b>44</b> |
|          | <b>ANNEX A: Flow diagram of the reference method</b>   | <b>46</b> |
|          | <b>ANNEX B: Flow diagram of the alternative method</b>   | <b>47</b> |
|          | <b>ANNEX C: Raw data from the sensitivity study</b>  | <b>48</b> |
|          | <b>ANNEX D: Raw data from the RLOD study</b>   | <b>65</b> |
|          | <b>ANNEX E: Raw data on inclusivity and exclusivity</b>  | <b>70</b> |

**ANNEX F: Collaborators in ILS**

**72**

## 1 Introduction

In this project a MicroVal validation study, based on ISO 16140-2:2016, of an alternative method for the detection of for the detection of Total Viable Count at a threshold of 10 cfu per g in a broad range of foods was carried out by Campden BRI as the MicroVal Expert Laboratory.

The alternative method used was:

The Total Viable Count (NF-TVC) vial (9 mL) is suitable for dilute-to-specification monitoring in a variety of food matrices. As organisms grow in the broth medium, the carbon dioxide (CO<sub>2</sub>) produced diffuses through a membrane layer into a soft agar plug containing a dye indicator. The color change in the dye is read by the Soleris® instrument.

The time to growth detection in the Soleris® system is correlated to the level of target organisms present in the sample, with higher levels of contamination having a shorter detection time. In this validation, the method was used to determine a defined threshold of product contamination of greater than 10 cfu per g. If successful, the alternative method would be used as a rapid screen for higher throughput product monitoring of a broad range of food products, reducing hands on analysis time as well as time to result.

The reference method used was:

ISO 4833-1:2013; Microbiology of the food chain -- Horizontal method for the enumeration of microorganisms -- Part 1: Colony count at 30°C by the pour plate technique

Although the reference method is generally used to enumerate the level of total viable count, in this validation it was used to establish if levels of the target organism exceeded the defined detection threshold of 10 CFU per g. Thus, the best approach to this semi-quantitative validation was to use a Qualitative presence/absence approach, but to set the presence/absence limit at 10 CFU/g.

Scope of the validation study is: a broad range of foods

Categories included:

- Dairy (heat processed)
- Fishery products (raw and RTC)
- Fresh and processed produce
- Multicomponent foods or meal components
- Meat and poultry products (Raw and RTC)

Criteria evaluated during the study have been:

- Method Comparison Study (MCS)
  - Sensitivity study
  - Relative level of detection study
  - Inclusivity and exclusivity study

The **method comparison study conclusions** are:

- The observed values for ND-PD for the individual categories and for all categories meet the acceptability limits (observed values  $\leq$  AL).
- The RLOD values (using the confirmed alternative method results) meet the acceptability limit, which is 2.5 for unpaired studies, for all categories tested.
- The alternative Neogen Soleris® NF-TVC detection method is selective and specific.

The **inter-laboratory study conclusions** are:

- The observed value for ND-PD is lower than the acceptability limits.
- There are no individual categories tested for the IL.
- The data and interpretations comply with the EN ISO 16140-2:2016 requirements.

This report corresponds to the method comparison study and gathers the observed data and interpretations according to the EN ISO 16140- 2:2016 standard and the MicroVal technical committee interpretations v2.5.



## **2 Method protocols**

The Method Comparison Study was carried out using 10 gram portions of sample material.

Although the Reference and Alternative methods were performed with the same test portion, the analysis does not share a common enrichment step, due to the requirement for pH adjustment in the alternative method. In addition, each method used in the study followed different principles. The Reference method detects growth of the target organisms in agar as part of a pour plate protocol. The Alternative method enables growth of the target organisms in a liquid medium to the set threshold of detection. Due to the differences in methodology, the validation was carried out as an unpaired study.

### **2.1 Reference method**

See the flow diagram in Annex A.

Sample preparations used in the reference method and the alternative method were done according to ISO 6887-series.

### **2.2 Alternative method**

See the flow diagram of the alternative method in Annex B.

The alternative method principle is based on optical detection of microbial growth through the use of a medium containing dye indicators. As organisms grow in the broth medium, the carbon dioxide (CO<sub>2</sub>) produced diffuses through a membrane layer into a soft agar plug containing a dye indicator. The color change in the dye is read by the Soleris® instrument.

At defined time intervals during incubation, the soft agar plug at the bottom of the vial is monitored for colour changes which indicated growth of the target organism. Readings were collated by the software within the instrument that automatically analysed the data to provide the final results. An incubation time of 30 hours was selected as the incubation time for the test, although the Soleris® system was run for 48 hours to check for any slow growing strains. All samples with a detection time of < 30 hours were reported as detected, and any detection times over 30 hours have been reported as not detected.

Following incubation of the vial for the required time, positive vials were streaked to PCA and incubated for 24-72 hours at 30°C±1°C to check for the presence of colonies. The shortest incubation time of 24h used, unless otherwise specified. Confirmation was performed for the purpose of this certification study to verify positive results obtained in the vials, as opposed to a requirement of customers where confirmation is optional. Customers can perform confirmation for troubleshooting purposes. Streaking the vial on PCA can allow customers to isolate colonies for further identification.

### **2.3 Study design**

The Method Comparison Study was carried out using 10 gram test portions of sample material.

Unless otherwise stated, sample preparations were done according to ISO 6887-series.

See Table 1 for specific preparations used in the validation study.

### 3 Method comparison study

#### 3.1 Sensitivity Study

*The sensitivity study (SE) is the ability of the method selected to detect the analyte by either the reference or the alternative method.*

##### 3.1.1 Categories and sample types

A total of 5 Categories were included in this validation study.

A minimum of 60 Items for each Category were tested by both the reference method and the alternative method in the sensitivity study, with a minimum of 30 positive samples per Category.

Each Category was made up of 3 Types, with at least 20 Items representative for that Type.

In the ideal situation, 10 samples (50 %) tested per type should be positive and 10 negative, but should range between 25 % and 75 %. For each category, at least 30 samples will have a positive result by the reference and/or the alternative method.

NF-TVC vial is a method for the detection of total viable count. The categories tested in the study were anticipated to be naturally contaminated with organisms capable of growing in a total viable count. During the study, a serial dilution series was prepared for each food item, and 3-4 dilutions within the dilution series were tested by both the reference and alternative methods. A single dilution (e.g. -6) was selected per category using pre-defined criteria to carry forward for statistical analysis. To enable compliance with the requirements of ISO 16140-2, a dilution that gave between 25%-75% positive samples for each type by the reference method for each sample type. The approach used in the current validation followed the study design for a previously certified TVC vial test, study reference 2021LR94.

On the occasion where two dilutions fulfil this criteria, the plate with the greater dilution factor was selected for the following reasons:

- To ensure that the correct number of positive and negative samples are achieved
- Higher dilutions have a lower count per plate (typically 1-10 cfu), which is closer to the limit of detection and is therefore a greater challenge for the method

The categories, the types and the number of samples analyzed are presented in Table 1.

**Table 1 - Categories, types and number of samples analyzed**

| Category                            | Types                      | Items (examples)                     | Preparations | Number of samples analysed |
|-------------------------------------|----------------------------|--------------------------------------|--------------|----------------------------|
| Pasteurised milk and dairy products | Pasteurised milk           | Pasteurised milk                     | 6887-5       | 20                         |
|                                     | Pasteurised dairy products | Ice-cream, milk based drinks, cheese | 6887-5       | 20                         |

| Category   | Types  | Items (examples)  | Preparations   | Number of samples analysed |
|--|--|---|----------------|----------------------------|
|  | Dry milk products                                    | Milk powder, dessert powder   | 6887-1, 6887-4 | 20                         |
| Raw Fishery products                                       | Raw fish (unprocessed)                               | Raw salmon filet, tuna  | 6887-3         | 25                         |
|  | Ready to cook  | Fishcakes   | 6887-3         | 25                         |
|  | Crustaceans  | Shrimp, crab  | 6887-3         | 20                         |
| Produce and fruits (combined category fresh and processed) | Cut ready-to-eat vegetables/leafy greens and sprouts | Bagged pre-cut lettuce shredded carrot, radish sprouts, alfalfa       | 6887-4         | 20                         |
|  | Fresh fruit/Cut RTE fruit                            | Cut fruits, freshly squeezed juice, smoothies                         | 6887-4         | 20                         |
|  | Heat treated fruit and vegetables                    | Pasteurised smoothies/juice, blanched frozen vegetables               | 6887-4         | 20                         |
| Multi-component foods or meal components                   | Composite foods with substantial raw ingredients     | Chilled pasta salad, sandwiches                                       | 6887-2         | 20                         |
|  | RTRH/RTE foods (chilled, frozen)                     | Cooked chilled pasta, frozen fries, rice products,                    | 6887-2         | 20                         |
|  | Mayonnaise based deli-salads                         | Vegetable salad   | 6887-2         | 20                         |
| Raw and Ready to cook RTC Meat and poultry                 | Raw poultry and meat cuts                            | Raw chicken, beef, pork, turkey                                       | 6887-2         | 20                         |
|  | Raw processed meat                                   | Frozen burger patties, pork meat balls, seasoned raw meat, lamb mince | 6887-2         | 20                         |
|  | RTC processed poultry                                | seasoned chicken, turkey meat balls                                   | 6887-2         | 20                         |

310 samples were analyzed. The distribution of positive and negative samples per tested category and type is given respectively in Table 2.

**Table 2 - Distribution per tested category and type**

| Category    |  | Type  |  | Number of positives* | Number of negatives |
|-------------|--|-------|--|----------------------|---------------------|
| 1           | Milk and dairy products                    | a     | Pasteurised milk                                     | 7                    | 13                  |
|             |  | b     | Pasteurised dairy products                           | 9                    | 11                  |
|             |  | c     | Dry milk products                                    | 14                   | 6                   |
|             |  | Total |  | 30                   | 30                  |
| 2           | Fishery products (raw and RTC)             | a     | Raw fish (unprocessed)                               | 16                   | 9                   |
|             |  | b     | RTC fish   | 12                   | 13                  |
|             |  | c     | Crustaceans  | 9                    | 11                  |
|             |  | Total |  | 37                   | 33                  |
| 3           | Fresh produce and fruits                   | a     | Cut ready-to-eat vegetables/leafy greens and sprouts | 13                   | 7                   |
|             |  | b     | Fresh fruit/cut RTE fruit and vegetable products     | 14                   | 6                   |
|             |  | c     | Heat treated fruit and vegetables                    | 8                    | 12                  |
|             |  | Total |  | 35                   | 25                  |
| 4           | Multi-component foods or meal components   | a     | Composite foods with substantial raw ingredients     | 11                   | 9                   |
|             |  | b     | RTRH/RTE foods (chilled, frozen)                     | 13                   | 7                   |
|             |  | c     | Mayonnaise based deli salads                         | 7                    | 13                  |
|             |  | Total |  | 31                   | 29                  |
| 5           | Raw and ready to cook RTC meat and poultry | a     | Raw poultry and meat cuts                            | 9                    | 11                  |
|             |  | b     | Raw processed meat                                   | 11                   | 9                   |
|             |  | c     | RTC processed meat and poultry                       | 14                   | 6                   |
|             |  | Total |  | 34                   | 26                  |
| Grand total |  |       |  | 167                  | 143                 |

\*Positive by at least one of the methods

### 3.1.2 Test sample preparation

100 % of the samples analysed in the study were naturally contaminated therefore artificial contaminations were not needed.

### 3.1.3 Confirmation protocols

Following 48 hours of incubation in the Soleris® instrument, each vial was streaked to PCA and incubated at 30±1°C for 24 hours to check for visible colonies.

### 3.1.4 Sensitivity study results

All raw data on the sensitivity study are given in Annex C.

Table 3 shows the summary of results of the reference method and the alternative methods for **all Categories**.

Table 4 shows the Interpretation of sample results between the reference and alternative method (based on the confirmed alternative method).

**Table 3 - Summary of sensitivity study results – all categories**

|   | Reference method<br>positive (R+)             | Reference method<br>negative (R-)             |
|---|---|---|
| <b>Alternative method<br/>positive (A+)</b> | Positive agreement (R+/A+)<br><b>PA = 134</b> | Positive deviation (R-/A+)<br><b>PD = 17</b>  |
| <b>Alternative method<br/>negative (A-)</b> | Negative deviation (R+/A-)<br><b>ND = 18</b>  | Negative agreement (R-/A-)<br><b>NA = 143</b> |

**Table 4 – Interpretation of sample results between the reference and alternative method (based on the confirmed alternative method)**

| Category |  | Type  |  | PA | NA <sup>1</sup> | PD | ND <sup>2</sup> | PPNA <sup>3</sup> | PPND <sup>3</sup> |
|----------|--|-------|--|----|-----------------|----|-----------------|-------------------|-------------------|
| 1        | Milk and dairy products                  | a     | Pasteurised milk                                     | 7  | 13              | 0  | 0               | 0                 | 0                 |
|          |  | b     | Pasteurised dairy products                           | 8  | 11              | 1  | 0               | 0                 | 0                 |
|          |  | c     | Dry milk products                                    | 14 | 6               | 0  | 0               | 0                 | 0                 |
|          |  | Total |  | 29 | 30              | 1  | 0               | 0                 | 0                 |
| 2        | Fishery products (raw and RTC)           | a     | Raw fish (unprocessed)                               | 12 | 9               | 1  | 3               | 1                 | 0                 |
|          |  | b     | RTC fish   | 10 | 13              | 0  | 2               | 0                 | 0                 |
|          |  | c     | Crustaceans  | 4  | 11              | 3  | 2               | 0                 | 0                 |
|          |  | Total |  | 26 | 33              | 4  | 7               | 1                 | 0                 |
| 3        | Fresh produce and fruits                 | a     | Cut ready-to-eat vegetables/leafy greens and sprouts | 12 | 7               | 1  | 0               | 0                 | 0                 |
|          |  | b     | Fresh fruit/cut RTE fruit and vegetable products     | 13 | 6               | 1  | 0               | 0                 | 0                 |
|          |  | c     | Heat treated fruit and vegetables                    | 2  | 12              | 1  | 5               | 0                 | 0                 |
|          |  | Total |  | 27 | 25              | 3  | 5               | 0                 | 0                 |
| 4        | Multi-component foods or meal components | a     | Composite foods with substantial raw ingredients     | 6  | 9               | 2  | 3               | 0                 | 0                 |
|          |  | b     | RTRH/RTE foods (chilled, frozen)                     | 13 | 7               | 0  | 0               | 0                 | 0                 |
|          |  | c     | Mayonnaise based deli salads                         | 6  | 13              | 0  | 1               | 0                 | 0                 |
|          |  | Total |  | 25 | 29              | 2  | 4               | 0                 | 0                 |

| Category    |  | Type  |                                | PA  | NA <sup>1</sup> | PD | ND <sup>2</sup> | PPNA <sup>3</sup> | PPND <sup>3</sup> |
|-------------|--|-------|--------------------------------|-----|-----------------|----|-----------------|-------------------|-------------------|
| 5           | Raw and ready to cook RTC meat and poultry | a     | Raw poultry and meat cuts      | 7   | 11              | 1  | 1               | 0                 | 0                 |
|             |  | b     | Raw processed meat             | 7   | 9               | 3  | 1               | 0                 | 0                 |
|             |  | c     | RTC processed meat and poultry | 13  | 6               | 1  | 0               | 0                 | 0                 |
|             |  | Total |                                | 27  | 26              | 5  | 2               | 0                 | 0                 |
| Grand total |  |       |                                | 134 | 143             | 15 | 18              | 1                 | 0                 |

<sup>1</sup> NA: including PPNA, <sup>2</sup> ND: including PPND, <sup>3</sup> FP = PPNA + PPND

### 3.1.5 Sensitivity study calculations

The sensitivity study parameters as specified in Table 5 were calculated for all Categories and Types, and the overview is given in Table 6.

**Table 5 – Formula to calculate the sensitivity parameters**

|   |  |
|---|--|
| Sensitivity for the alternative method          | $SE_{alt} = \frac{(PA + PD)}{(PA + ND + PD)} \times 100\%$ |
| Sensitivity for the reference method            | $SE_{ref} = \frac{(PA + ND)}{(PA + ND + PD)} \times 100\%$ |
| Relative trueness                               | $RT = \frac{(PA + NA)}{N} \times 100\%$                    |
| False positive ratio for the alternative method | $FPR = \frac{(FP)}{NA} \times 100\%$                       |



**Table 6 - Overview calculated sensitivity parameters per Category and Type**

| Category    |  | Type  |  | PA  | NA  | PD | ND | PPNA | PPND | Total | SE alt | SE ref | RT  | FPR |
|-------------|--|-------|--|-----|-----|----|----|------|------|-------|--------|--------|-----|-----|
| 1           | Milk and dairy products                    | a     | Pasteurised milk                                     | 7   | 13  | 0  | 0  | 0    | 0    | 20    | 100    | 100    | 100 | 0   |
|             |  | b     | Pasteurised dairy products                           | 8   | 11  | 1  | 0  | 0    | 0    | 20    | 100    | 100    | 100 | 0   |
|             |  | c     | Dry milk products                                    | 14  | 6   | 0  | 0  | 0    | 0    | 20    | 100    | 100    | 100 | 0   |
|             |  | Total |  | 29  | 30  | 1  | 0  | 0    | 0    | 60    | 100    | 100    | 100 | 0   |
| 2           | Fishery products (raw and RTC)             | a     | Raw fish (unprocessed)                               | 12  | 9   | 1  | 3  | 1    | 0    | 25    | 92     | 75     | 69  | 11  |
|             |  | b     | RTC fish   | 10  | 13  | 0  | 2  | 0    | 0    | 25    | 100    | 87     | 87  | 0   |
|             |  | c     | Crustaceans  | 4   | 11  | 3  | 2  | 0    | 0    | 20    | 100    | 85     | 88  | 0   |
|             |  | Total |  | 26  | 33  | 4  | 7  | 1    | 0    | 70    | 98     | 83     | 82  | 3   |
| 3           | Fresh produce and fruits                   | a     | Cut ready-to-eat vegetables/leafy greens and sprouts | 12  | 7   | 1  | 0  | 0    | 0    | 20    | 100    | 100    | 100 | 0   |
|             |  | b     | Fresh fruit/cut RTE fruit and vegetable products     | 13  | 6   | 1  | 0  | 0    | 0    | 20    | 100    | 100    | 100 | 0   |
|             |  | c     | Heat treated fruit and vegetables                    | 2   | 12  | 1  | 5  | 0    | 0    | 20    | 100    | 71     | 72  | 0   |
|             |  | Total |  | 27  | 25  | 3  | 5  | 0    | 0    | 60    | 100    | 83     | 85  | 0   |
| 4           | Multi-component foods or meal components   | a     | Composite foods with substantial raw ingredients     | 6   | 9   | 2  | 3  | 0    | 0    | 20    | 100    | 75     | 79  | 0   |
|             |  | b     | RTRH/RTE foods (chilled, frozen)                     | 13  | 7   | 0  | 0  | 0    | 0    | 20    | 100    | 100    | 100 | 0   |
|             |  | c     | Mayonnaise based deli salads                         | 6   | 13  | 0  | 1  | 0    | 0    | 20    | 100    | 93     | 93  | 0   |
|             |  | Total |  | 25  | 29  | 2  | 4  | 0    | 0    | 60    | 100    | 88     | 89  | 0   |
| 5           | Raw and ready to cook RTC meat and poultry | a     | Raw poultry and meat cuts                            | 6   | 11  | 2  | 1  | 0    | 0    | 20    | 100    | 92     | 93  | 0   |
|             |  | b     | Raw processed meat                                   | 7   | 9   | 3  | 1  | 0    | 0    | 20    | 100    | 90     | 92  | 0   |
|             |  | c     | RTC processed meat and poultry                       | 13  | 6   | 1  | 0  | 0    | 0    | 20    | 100    | 100    | 100 | 0   |
|             |  | Total |  | 26  | 26  | 6  | 2  | 0    | 0    | 60    | 100    | 93     | 94  | 0   |
| Grand total |  |       |  | 133 | 143 | 16 | 18 | 1    | 0    | 310   | 99     | 89     | 90  | 0.7 |

<sup>1</sup> NA: including PPNA, <sup>2</sup> ND: including PPND, <sup>3</sup>FP = PPNA + PPND



### 3.1.6 Discordant results

Negative deviations are listed in Table 7.

**Table 7 - Negative deviations**

| Category                                   | Type   | Sample                                   | Sample no | Plate count for reference method |
|--|--|--|-----------|----------------------------------|
| Fishery products (raw and RTC)             | Raw fish (unprocessed)                           | Scottish tail salmon fillets             | S42       | 1                                |
| Fishery products (raw and RTC)             | Raw fish (unprocessed)                           | Sea bass fillets                         | T17       | 3                                |
| Fishery products (raw and RTC)             | Raw fish (unprocessed)                           | Raw squid rings                          | F62       | 9                                |
| Fishery products (raw and RTC)             | RTC fish   | Creamy cod fishcakes                     | F7        | 2                                |
| Fishery products (raw and RTC)             | RTC fish   | Smoked haddock fishcakes                 | F13       | 6                                |
| Fishery products (raw and RTC)             | Crustaceans                                      | White crab 2                             | S5        | 1                                |
| Fishery products (raw and RTC)             | Crustaceans                                      | White & Brown Crab                       | S6        | 1                                |
| Fresh produce and fruits                   | Heat treated fruit and vegetables                | Orange juice                             | T42       | 1                                |
| Fresh produce and fruits                   | Heat treated fruit and vegetables                | Mango baby food                          | P46       | 1                                |
| Fresh produce and fruits                   | Heat treated fruit and vegetables                | Apple, mango and banana baby food        | P47       | 1                                |
| Fresh produce and fruits                   | Heat treated fruit and vegetables                | Pear baby food                           | P49       | 1                                |
| Fresh produce and fruits                   | Heat treated fruit and vegetables                | V8 vegetable juice                       | P60       | 2                                |
| Multi-component foods or meal components   | Composite foods with substantial raw ingredients | Cheese and onion sandwich                | T50       | 2                                |
| Multi-component foods or meal components   | Composite foods with substantial raw ingredients | Tomato and basil pasta salad             | C18       | 4                                |
| Multi-component foods or meal components   | Composite foods with substantial raw ingredients | Honey and mustard chicken pasta          | C19       | 5                                |
| Multi-component foods or meal components   | Mayonnaise based deli salads                     | Chicken caesar salad                     | T58       | 2                                |
| Raw and ready to cook RTC meat and poultry | Raw poultry and meat cuts                        | Fresh British turkey breast mini fillets | R10       | 1                                |
| Raw and ready to cook RTC meat and poultry | Raw processed meat                               | Pork shoulder in bbq sauce               | T7        | 2                                |

Positive deviations are listed in Table 8.





**Table 8 - Positive deviations**

| Category                                   | Type   | Sample                                | Sample n° |
|--|--|---------------------------------------|-----------|
| Fishery products (raw and RTC)             | Raw fish (unprocessed)                               | Salmon fillet                         | S41       |
| Fishery products (raw and RTC)             | Crustaceans  | Cold water large prawns               | T26       |
| Fishery products (raw and RTC)             | Crustaceans  | White crab                            | S1        |
| Fishery products (raw and RTC)             | Crustaceans  | Premium crab meat                     | S4        |
| Fresh produce and fruits                   | Cut ready-to-eat vegetables/leafy greens and sprouts | Rocket salad                          | P20       |
| Fresh produce and fruits                   | Fresh fruit/cut RTE fruit and vegetable products     | Mango                                 | P29       |
| Fresh produce and fruits                   | Fresh fruit/cut RTE fruit and vegetable products     | Tropical juice                        | P57       |
| Multicomponent foods or meal components    | Composite foods with substantial raw ingredients     | Cheese and pickle sandwich            | C13       |
| Multi-component foods or meal components   | Composite foods with substantial raw ingredients     | tesco chicken and chorizo pasta salad | C14       |
| Raw and ready to cook RTC meat and poultry | Raw poultry and meat cuts                            | Turkey thigh mince 7% fat             | T4        |
| Raw and ready to cook RTC meat and poultry | Raw processed meat                                   | Fire pit sweet and smoky beef kebabs  | T8        |
| Raw and ready to cook RTC meat and poultry | Raw processed meat                                   | Quarter pounder burgers               | R37       |
| Raw and ready to cook RTC meat and poultry | Raw processed meat                                   | Venison burgers                       | R40       |
| Raw and ready to cook RTC meat and poultry | RTC processed meat and poultry                       | British beef quarter pounders         | R58       |
| Milk and dairy products                    | Pasteurised dairy products                           | Soured cream                          | T67       |

There were 18 negative deviations, and 15 positive deviations observed in the sensitivity study. Further analysis revealed that the negative deviations were seen in 4 out of the 5 categories tested. Statistical analysis was performed to check the potential possibility of a link between the dilution used and the deviations obtained. The results demonstrated that there was no significant difference between the average dilutions available between negative and positive deviations ( $P=0.4$ ). For all negative deviations, the level obtained on the reference method is between 1-9 cfu/plate. The method shows equivalent performance to the reference method at low levels in the RLOD portion of the study.

The ND-PD values were highest in the multicomponent foods and produce and fruits categories. Similar levels of negative and positive deviations were noted in the other three categories tested.

Most of the negative deviations seen in the produce and fruits category were obtained with heat processed fruit and vegetable samples. Vegetable juice was also tested in the RLOD study and achieved an RLOD of 0.640. The RLOD value for the vegetable juice was below the AL of 2.5 for an unpaired study and does not indicate systematic bias with the sample type. One possible explanation for the deviations could be the presence of moulds, as mould colonies were observed in several heat processed fruit and vegetable samples tested. All 5 yeast and mould strains analysed in the inclusivity study were detected within 30 hours, although the incubation time and temperature is not optimised for yeasts and moulds. It is recommended to use a test specific to yeast and moulds if only yeast and moulds are expected in the samples to be analysed.

There were 11 samples (10 PA, 1 PD) that required the incubation time of the confirmation plate to be increased from 24 hours to 72 hours. These samples are identified in the raw data, and are from the fishery product category. One possible explanation for this is the organisms present in these samples are slower growing.

The analysis of discordant results according to ISO 16140-2:2016 for an unpaired study is given in Table 9.

**Table 9 – Interpretation of the sensitivity study results (unpaired study)**

| Category   | Negative Deviations (ND <sup>1</sup> ) | Positive Deviations (PD) | ND-PD    | Acceptability Limit (AL) |
|--|--|--------------------------|----------|--------------------------|
| Milk and dairy products                                | 0                                      | 1                        | -1       | 3                        |
| Fishery products (raw and RTC)                         | 7                                      | 4                        | 3        | 3                        |
| Produce and fruits (combined fresh and heat processed) | 5                                      | 3                        | 2        | 3                        |
| Multi-component foods or meal components               | 4                                      | 2                        | 2        | 3                        |
| Raw and ready to cook RTC meat and poultry             | 2                                      | 5                        | -3       | 3                        |
| <b>Total</b>   | <b>18</b>                              | <b>15</b>                | <b>3</b> | <b>5</b>                 |

<sup>1</sup> ND: including PPND

### 3.1.7 Conclusion sensitivity study

The observed values for ND-PD for the individual categories and for all categories meet the acceptability limits (observed values  $\leq$  AL).

## 3.2 Relative level of detection study

*The relative level of detection is the level of detection at  $P = 0,50$  ( $LOD_{50}$ ) of the alternative method divided by the level of detection at  $P = 0,50$  ( $LOD_{50}$ ) of the reference method.*

### 3.2.1 Categories, sample types and strains

One sample type and one relevant target micro-organism for this sample type was chosen for each of the Categories in this validation study, as shown in Table 10. Canned products were selected for the RLOD study to enable a consistent low level of artificial contamination to be achieved.

**Table 10 - List of selected types and strains per category, as tested within the relative level of detection study.**

| Category   | Type  | Item                         | Strain   | Seeding or spiking procedure   | Log injury |
|--|---|------------------------------|--|--|------------|
| Dairy products (heat processed)                            | Dry dairy products                              | Milk powder                  | <i>Bacillus cereus</i> CRA 1724 isolated from Dried milk         | Seeding – inoculated with a lyophilized culture and stored for >2 weeks at ambient | N/A        |
| Fishery products (raw and RTC)                             | Canned ambient stable fish                      | Tuna Chunks in Sunflower Oil | <i>Pseudomonas fragi</i> CRA7222 Isolated from spoiled fish      | Spiking – culture heat stressed at 55°C for 15 minutes                             | 0.57       |
| Produce and fruits (combined category fresh and processed) | Heat processed fruit and vegetables             | Vegetable juice              | <i>E.coli</i> CRA3379 isolated from spinach                      | Spiking – culture heat stressed at 55°C for 15 minutes                             | 0.28       |
| Raw and RTC Meat and poultry (Combined category)           | Canned ambient stable meat                      | Chopped Ham & Pork           | <i>Staphylococcus warneri</i> CRA 3198 Isolated from dry sausage | Spiking – culture heat stressed at 55°C for 15 minutes                             | 0.77       |
| Multicomponent foods or meal components                    | Ready to (re)heat food: ambient stable (canned) | Spaghetti hoops              | <i>Hafnia alvei</i> CRA7417 isolated from pate                   | Spiking – culture heat stressed at 55°C for 15 minutes                             | 0.65       |

### 3.2.2 Test sample preparations

Three levels of artificial contamination were prepared for each type:

- Negative control level: One non-inoculated in order to get 5 test portions,
- Low level: One inoculated between 8 and 12 CFU/g in order to get 20 test portions providing fractional recovery,
- Higher level: One inoculated between 15 and 25 CFU/sample in order to get 5 test portions contaminated at a higher level.

### 3.2.3 RLOD study results

The tabulated raw data on the RLOD study are given in Annex D.

The RLOD calculations were performed using the Excel spread sheet (version 3 15-08-2015) of the international standard as described in ISO 16140-2: 2016.

The RLOD per Category is given in Table 11.

**Table 11 – Presentation of RLOD before and after confirmation of the alternative method results**

| Type (Category)  | RLOD using the alternative method results | RLOD using the confirmed alternative method results |
|--|---|---|
| Dairy products (heat processed)                            | 0.838                                     | 0.838   |
| Fishery products (raw and RTC)                             | 1.000                                     | 1.513   |
| Produce and fruits (combined category fresh and processed) | 0.640                                     | 0.640   |
| Raw and RTC Meat and poultry                               | 0.456                                     | 0.456   |
| Multicomponent foods or meal components                    | 1.161                                     | 1.161   |
| <b>Combined</b>  | <b>0.793</b>                              | <b>0.854</b>  |

The RLOD values (using the confirmed alternative method results) meet the acceptability limit, 2.5 for unpaired studies, for all categories tested.

In addition, LOD50 values were calculated using the equations quoted in Wilrich and Wilrich (2009) Journal of AOAC International 92 (6) 1763-1772 downloaded from

[www.wiwiw.fu-berlin.de/fachbereich/vwl/iso/ehemalige/wilrich.index.html](http://www.wiwiw.fu-berlin.de/fachbereich/vwl/iso/ehemalige/wilrich.index.html)

The LOD50 per Category is given in Tables 12 and 13.

**Table 12 – LOD<sub>50</sub> calculation for presumptive results**

| Results of the PODLOD calculations - with the LOD relating to <i>d</i> in cfu / g |  |                        |                            |                                     |  |                                   |                                   |  |                                    |                                    |   |
|---|--|------------------------|----------------------------|-------------------------------------|--|-----------------------------------|-----------------------------------|--|------------------------------------|------------------------------------|---|
| Matrix  |  | Matrix effect<br>$F_i$ | Log matrix effect<br>$f_i$ | SD of log matrix effect<br>$s_{fi}$ | LOD <sub>50%</sub> = 50% limit of detection in cfu / g |                                   |                                   | LOD <sub>95%</sub> = 95% limit of detection in cfu / g |                                    |                                    | Test statistic matrix effect<br>$ z_i $ |
| No. <i>i</i>  | Designation Matrix <sub><i>i</i></sub> |                        |                            |                                     | Detection limit<br>$d_{0.5,i}$                         | Lower conf. limit<br>$d_{0.5,iL}$ | Upper conf. limit<br>$d_{0.5,iU}$ | Detection limit<br>$d_{0.95,i}$                        | Lower conf. limit<br>$d_{0.95,iL}$ | Upper conf. limit<br>$d_{0.95,iU}$ |   |
| 1   | Dairy                                  | 0.010                  | -4.567                     | 0.259                               | 6.671  | 3.977                             | 11.191                            | 28.831   | 17.186                             | 48.366                             | 2.9E-21                                 |
| 2   | Fishery products                       | 0.012                  | -4.399                     | 0.256                               | 5.638  | 3.382                             | 9.400                             | 24.367   | 14.616                             | 40.625                             | 2.1E-23                                 |
| 3   | Vegetable juice                        | 0.013                  | -4.317                     | 0.257                               | 5.194  | 3.110                             | 8.676                             | 22.449   | 13.439                             | 37.497                             | 2.8E-21                                 |
| 4   | Meat and poultry                       | 0.008                  | -4.772                     | 0.263                               | 8.192  | 4.842                             | 13.860                            | 35.405   | 20.926                             | 59.901                             | 3.1E-21                                 |
| 5   | Multicomponent foods                   | 0.021                  | -3.860                     | 0.286                               | 3.289  | 1.856                             | 5.828                             | 14.213   | 8.020                              | 25.190                             | 2.5E-21                                 |
| Combined results  |  | 0.012                  | -4.413                     | 0.115                               | 5.717  | 4.543                             | 7.194                             | 24.708   | 19.636                             | 31.091                             | 5.6E-21                                 |
| The combined results are based on the data of matrices 1, 2, 3, 4 and 5           |  |                        |                            |                                     |  |                                   |                                   |  |                                    |                                    |   |

**Table 13 - LOD<sub>50</sub> calculation for confirmed results**

| Results of the PODLOD calculations - with the LOD relating to <i>d</i> in cfu / g |  |                        |                            |                                     |  |                                   |                                   |  |                                    |                                    |   |
|---|--|------------------------|----------------------------|-------------------------------------|--|-----------------------------------|-----------------------------------|--|------------------------------------|------------------------------------|---|
| Matrix  |  | Matrix effect<br>$F_i$ | Log matrix effect<br>$f_i$ | SD of log matrix effect<br>$s_{fi}$ | LOD <sub>50%</sub> = 50% limit of detection in cfu / g |                                   |                                   | LOD <sub>95%</sub> = 95% limit of detection in cfu / g |                                    |                                    | Test statistic matrix effect<br>$ z_i $ |
| No. <i>i</i>  | Designation Matrix <sub><i>i</i></sub> |                        |                            |                                     | Detection limit<br>$d_{0.5,i}$                         | Lower conf. limit<br>$d_{0.5,iL}$ | Upper conf. limit<br>$d_{0.5,iU}$ | Detection limit<br>$d_{0.95,i}$                        | Lower conf. limit<br>$d_{0.95,iL}$ | Upper conf. limit<br>$d_{0.95,iU}$ |   |
| 1   | Dairy                                  | 0.010                  | -4.567                     | 0.259                               | 6.671  | 3.977                             | 11.191                            | 28.831   | 17.186                             | 48.366                             | 2.9E-21                                 |
| 2   | Fishery products                       | 0.012                  | -4.399                     | 0.256                               | 5.638  | 3.382                             | 9.400                             | 24.367   | 14.616                             | 40.625                             | 2.1E-23                                 |
| 3   | Vegetable juice                        | 0.009                  | -4.668                     | 0.260                               | 7.384  | 4.393                             | 12.411                            | 31.912   | 18.986                             | 53.639                             | 3.0E-21                                 |
| 4   | Meat and poultry                       | 0.008                  | -4.772                     | 0.263                               | 8.192  | 4.842                             | 13.860                            | 35.405   | 20.926                             | 59.901                             | 3.1E-21                                 |
| 5   | Multicomponent foods                   | 0.021                  | -3.860                     | 0.286                               | 3.289  | 1.856                             | 5.828                             | 14.213   | 8.020                              | 25.190                             | 2.5E-21                                 |
| Combined results  |  | 0.011                  | -4.484                     | 0.115                               | 6.142  | 4.881                             | 7.728                             | 26.544   | 21.094                             | 33.402                             | 5.7E-21                                 |
| The combined results are based on the data of matrices 1, 2, 3, 4 and 5           |  |                        |                            |                                     |  |                                   |                                   |  |                                    |                                    |   |

### 3.2.4 Conclusion RLOD study

The RLOD values (using the confirmed alternative method results) meet the acceptability limit, 2.5 for unpaired studies, for all categories tested.

### 3.3 Inclusivity/exclusivity study

*Inclusivity is the ability of the alternative method to detect the target analyte from a wide range of strains.*

*Exclusivity is the lack of interference from a relevant range of non-target strains of the alternative method.*

#### 3.3.1 Protocols

**Inclusivity:** 50 strains of target microorganisms (bacteria, yeasts and moulds) able to grow in aerobic conditions at 30°C were freshly cultured in an appropriate non-selective broth at 30°C. Dilutions were made in order for the vials to be inoculated at 10-100 times greater than the minimal level of detection (10-100 cfu in vial). The alternative method protocol was then performed (no sample material was added).

No exclusivity study was carried out for this method, as the vial detects all organisms capable of growing in an aerobic plate count, hence there are no non-target organisms.

### *3.3.2 Results inclusivity and exclusivity study*

All raw data on inclusivity and exclusivity are given in Annex E.

A total of 50 strains were tested for **inclusivity**. All 50 isolates showed the expected positive result with detection times ranging from 10.9 to 29.6 hours.

### *3.3.3 Conclusion inclusivity and exclusivity study*

The alternative Neogen Soleris® NF-TVC detection method is selective and specific.

## **4 Conclusions Method Comparison Study**

Overall, the conclusions for the Method Comparison Study are:

The observed values for ND-PD for the individual categories and for all categories meet the acceptability limits (observed values  $\leq$  AL).

The RLOD values (using the confirmed alternative method results) meet the acceptability limit, which is 2.5 for unpaired studies, for all categories tested.

The alternative Neogen Soleris® NF-TVC detection method is selective and specific.

## 5 Interlaboratory study

The interlaboratory study is a study performed by multiple laboratories testing identical samples at the same time, the results of which are used to estimate alternative-method performance parameters.

### 5.1 Study organisation

#### *Collaborators number*

Samples were sent to 14 participants from 9 laboratories. The participants were located in seven countries: Switzerland, Spain, Netherlands, Ireland, France, United States and England. The participant details are shown in Annex F.

In addition to the 14 external participants, an internal participant from the analytical services team at Campden BRI, also took part in the study. This participant is in addition to the expert laboratory participant who was unfamiliar with the alternative method.

#### *Matrix and strain used*

A UHT soya-based dessert was inoculated with *Staphylococcus saprophyticus* CRA 314 isolated from a factory environment.

#### *Samples*

For each participant, 24 blind coded samples were prepared for detection of total viable count analysis by the reference and alternative methods.

Samples were prepared and inoculated on Monday 10<sup>th</sup> June 2024.

#### *Inoculation*

The targeted inoculation levels were the following:

- Level 0: 0 CFU/ g,
- Level 1: 11 CFU/ g, inoculation level providing as much as possible fractional positive recovery data;
- Level 2: 30 CFU/g

Following inoculation, samples were stored frozen for 24h before shipment to the participants.

#### *Labelling and shipping*

Blind coded samples were placed in isothermal boxes, which contained cooling blocks, and express-shipped to the different laboratories taking part in the study.

A frozen water vial was added to the package in order for participants to record the temperature on arrival.

A temperature probe sensor was added to the package in order to register the temperature profile during the transport, the package delivery and storage until analyses.

Samples were shipped in 24 h to 48h to the involved laboratories. The temperature conditions had to stay lower or equal to 1°C during transport and sample storage in the labs.

### Analyses

Collaborative study laboratories and the expert laboratory carried out the analyses on 17<sup>th</sup> June 2024 with the alternative and reference methods. The analyses by the reference method and the alternative method were performed on the same day.

## 5.2 Experimental parameters controls

### 5.2.1 Detection of total viable count in the matrix before inoculation

To detect the presence of total viable count the reference method was performed on eight portions (10 g) before the inoculation. All the results were negative indicating no background contamination in the samples.

### 5.2.2 Strain stability during transport

To test the stability of the strain in the matrix. 500cfu/g was inoculated into samples and tested after 7 days of storage, at two storage protocols, that mimic the conditions of storage and transit.

Storage protocol 1 – 24h frozen, 48h in temperature-controlled storage box, 4 days frozen

Storage protocol 2 – 24h chilled, 48h in temperature-controlled storage box, 4 days chilled

**Table 14 – Total viable count stability in the matrix**

| Day                        | Reference method (cfu/g) |             |
|----------------------------|--------------------------|-------------|
|                            | Replicate 1              | Replicate 2 |
| Day 0                      | 260                      | 150         |
| Storage protocol 1 (Day 7) | 100                      | 40          |
| Storage protocol 2 (Day 7) | 70                       | 80          |

The results show that there was a slight die-off after storage for 7 days.



Samples were then inoculated at lower levels to determine the detection of both methods before and after storage, the results are shown in Table 15.

**Table 15 – Stability results**

| Day                        | Reference method |       |          |       |         |       | Alternative method |       |          |       |         |       |
|----------------------------|------------------|-------|----------|-------|---------|-------|--------------------|-------|----------|-------|---------|-------|
|                            | 0cfu/g           |       | 10 cfu/g |       | 25cfu/g |       | 0cfu/g             |       | 10 cfu/g |       | 25cfu/g |       |
|                            | Rep 1            | Rep 2 | Rep 1    | Rep 2 | Rep 1   | Rep 2 | Rep 1              | Rep 2 | Rep 1    | Rep 2 | Rep 1   | Rep 2 |
| Day 0                      | -                | -     | -        | +     | +       | +     | -                  | -     | +        | -     | +       | +     |
| Storage protocol 1 (day 7) | -                | -     | +        | -     | +       | +     | -                  | -     | -        | -     | +       | +     |
| Storage protocol 2 (day 7) | +                | +     | +        | +     | +       | -     | -                  | -     | +        | +     | +       | +     |

The results show fractional positives with both methods at day 7. Storage protocol 1 was used as it appears to reduce the risk of positive uninoculated samples.

### 5.2.3 Logistic conditions

The temperatures measured at reception by the collaborators, the temperatures registered by the thermo-probe, and the receipt dates are given in Table 16.

**Table 16 - Sample temperatures at receipt**

| Collaborator | Temperature measured by the probe (°C) | Temperature measured at receipt (°C) | Receipt date and time | State of the package and samples at the receipt | Analysis date |
|--------------|--|--------------------------------------|-----------------------|---|---------------|
| 1            | Sample did not arrive                  |                                      |                       |   |               |
| 2            | -0.1                                   | Water control damaged                | 13/6/2024 14:20       | Ok – except for water control                   | 17/6/2024     |

| Collaborator | Temperature measured by the probe (°C) | Temperature measured at receipt (°C) | Receipt date and time   | State of the package and samples at the receipt | Analysis date |
|--------------|--|--------------------------------------|-------------------------|---|---------------|
| 3A           | -0.1                                   | N/A – frozen water control           | 14/6/2024 14:00         | Ok  | 17/6/2024     |
| 3B           | 0.9                                    | N/A – frozen water control           | 14/6/2024 14:00         | Ok  | 17/6/2024     |
| 4            | -1.1                                   | N/A – frozen water control           | 13/6/2024<br>Not listed | Ok  | 24/6/2024     |
| 5A           | -4.6                                   | N/A – frozen water control           | 13/6/2024 14:00         | Ok  | 17/6/2024     |
| 5B           | -1.1                                   | N/A – frozen water control           | 13/6/2024 14:00         | Ok  | 17/6/2024     |
| 6A           | -3.0                                   | N/A – frozen water control           | 13/6/2024 14:00         | Ok  | 17/6/2024     |
| 6B           | -1.6                                   | N/A – frozen water control           | 13/6/2024 14:00         | Ok  | 17/6/2024     |
| 7            | -1.0                                   | 3°C                                  | 13/6/2024 10:00         | Ok  | 17/6/2024     |
| 8A           | -1.4                                   | N/A – frozen water control           | 13/6/2024 14:50         | Ok  | 17/6/2024     |
| 8B           | -1.4                                   | N/A – frozen water control           | 13/6/2024 14:50         | Ok  | 17/6/2024     |
| 9A           | Data not provided – test not completed |                                      |                         |   |               |
| 9B           | Data not provided – test not completed |                                      |                         |   |               |

Participant 1 did not receive the samples due to customs issues. Participant 9A + 9B informed us that the samples were stored ambient on arrival and did not complete the testing.

Participant 4 received the samples on time, and stored them correctly, however, the set-up day was 7 days later than the other participants. This is outside the scope of the stability trials.

For the remainder of the participants, no problem was encountered during the transport or at receipt. Temperatures during shipment and at receipt were all correct.

### 5.3 Calculation and summary of data

#### 5.3.1 MicroVal Expert laboratory results

The results obtained by the expert laboratory are given in Table 17.

**Table 17 – Results obtained by the expert lab.**

| Level | Reference method | Alternative method |
|-------|------------------|--------------------|
| L0    | 0/8              | 0/8                |
| L1    | 6/8              | 6/8                |
| L2    | 8/8              | 8/8                |

#### 5.3.2 Results obtained by the collaborative laboratories

- Total Viable Count *detection*

12 out of 15 participants returned data for the study. The results obtained by the individual collaborators in the inter-laboratory study are summarised in Table 18 (reference method) and Table 19 (alternative method).

**Table 18 - Positive results by the reference method (ALL the collaborators)**

| Participant | Contamination level |     |     |
|-------------|---------------------|-----|-----|
|             | L0                  | L1  | L2  |
| 2           | 0/8                 | 6/8 | 8/8 |
| 3A          | 2/8                 | 6/8 | 8/8 |
| 3B          | 1/8                 | 7/8 | 8/8 |
| 4           | 5/8                 | 2/8 | 4/8 |
| 5A          | 3/8                 | 6/8 | 8/8 |
| 5B          | 2/8                 | 5/8 | 8/8 |

| Participant  | Contamination level          |                              |                              |
|--------------|------------------------------|------------------------------|------------------------------|
|              | L0                           | L1                           | L2                           |
| 6A           | 0/8                          | 6/8                          | 8/8                          |
| 6B           | 0/8                          | 6/8                          | 8/8                          |
| 7            | 0/8                          | 6/8                          | 8/8                          |
| 8A           | 0/8                          | 6/8                          | 8/8                          |
| 8B           | 0/8                          | 6/8                          | 8/8                          |
| 10           | 0/8                          | 6/8                          | 8/8                          |
| <b>TOTAL</b> | <b>P<sub>0</sub> = 13/96</b> | <b>P<sub>1</sub> = 69/96</b> | <b>P<sub>2</sub> = 92/96</b> |

**Table 19 - Positive results (before and after confirmation) by the alternative methods (ALL the collaborators)**

| Participant | Contamination level        |                           |                            |                           |                            |                           |
|-------------|----------------------------|---------------------------|----------------------------|---------------------------|----------------------------|---------------------------|
|             | L0                         |                           | L1                         |                           | L2                         |                           |
|             | <i>Before confirmation</i> | <i>After confirmation</i> | <i>Before confirmation</i> | <i>After confirmation</i> | <i>Before confirmation</i> | <i>After confirmation</i> |
| 2           | 0/8                        | 0/8                       | 6/8                        | 6/8                       | 8/8                        | 8/8                       |
| 3A          | 0/8                        | 0/8                       | 6/8                        | 6/8                       | 8/8                        | 8/8                       |
| 3B          | 0/8                        | 0/8                       | 6/8                        | 6/8                       | 8/8                        | 8/8                       |
| 4           | 1/8                        | 1/8                       | 6/8                        | 6/8                       | 8/8                        | 8/8                       |
| 5A          | 0/8                        | 0/8                       | 6/8                        | 6/8                       | 8/8                        | 8/8                       |
| 5B          | 0/8                        | 0/8                       | 6/8                        | 6/8                       | 8/8                        | 8/8                       |
| 6A          | 0/8                        | 0/8                       | 6/8                        | 6/8                       | 8/8                        | 8/8                       |

| Participant  | Contamination level            |                               |                                |                               |                                |                               |
|--------------|--------------------------------|-------------------------------|--------------------------------|-------------------------------|--------------------------------|-------------------------------|
|              | L0                             |                               | L1                             |                               | L2                             |                               |
|              | <i>Before<br/>confirmation</i> | <i>After<br/>confirmation</i> | <i>Before<br/>confirmation</i> | <i>After<br/>confirmation</i> | <i>Before<br/>confirmation</i> | <i>After<br/>confirmation</i> |
| 6B           | 0/8                            | 0/8                           | 6/8                            | 6/8                           | 8/8                            | 8/8                           |
| 7            | 0/8                            | 2/8                           | 0/8                            | 6/8                           | 0/8                            | 8/8                           |
| 8A           | 0/8                            | 0/8                           | 6/8                            | 6/8                           | 8/8                            | 8/8                           |
| 8B           | 0/8                            | 0/8                           | 6/8                            | 6/8                           | 8/8                            | 8/8                           |
| 10           | 0/8                            | 0/8                           | 6/8                            | 6/8                           | 8/8                            | 8/8                           |
| <b>TOTAL</b> | <b>P<sub>0</sub> = 1/96</b>    | <b>CP<sub>0</sub> = 3/8</b>   | <b>P<sub>1</sub> = 66/96</b>   | <b>CP<sub>1</sub> = 72/96</b> | <b>P<sub>2</sub> = 88/96</b>   | <b>CP<sub>2</sub> = 96/96</b> |

Analysis of the data revealed issues were experienced in sample analysis for 2 of the participants (4 and 7) and these are detailed below.

#### Participant 4

It was noted that participant 4 stored their samples frozen for 2 weeks instead of 1 week stated in the ILS instructions. A 2-week storage time is beyond the scope of the stability trial which was carried out for 1 week. Data analysis showed that storage has potential impact on reference method. The number of positives for each method is shown below:

| Level of contamination | Number of positives reference method | Number of positives alternative method |
|------------------------|--------------------------------------|--|
| Level 1                | 2                                    | 6                                      |
| Level 2                | 4                                    | 8                                      |

There are 5 positive blank samples for the reference method and 1 positive blank sample for the alternative method. Further investigations on the identification of the positive blank samples were inconclusive.

Due to the results of the investigations and the potential impact on the reference method, the results for participant 4 have been excluded from the analysis.

### Participant 7

There was no detection time recorded with the Soleris instrument for any samples tested by participant 7. On visual inspection, the vials showed visible signs of growth (e.g. change of plug colour, turbidity in broth contained in the vial). In addition, the vials were streaked and fractional positive confirmation results were received.

Further investigation found that the graphs were not plotting, which indicated that the instrument was not recording measurements. The instrument was out of calibration and using a previous software version.

Due to the instrument issues, the results for participant 7 were not taken forward for the analysis.

### Positive blank samples

For the remaining 10 participants in the study, four participants had up to 3 positive samples for the reference method. This translates to 10% positive results for the blank samples tested by the 10 participants taken forward for the analysis. Further investigations were carried out by the participants, including MALDI-ToF confirmations and/or images of the plates. Based on the results of these investigations, all remaining 10 participants have been included in the statistical analysis.

### Conclusion

Participants 4 and 7 have been excluded. The remaining labs with positive blank samples have been retained due to the results of the investigations. Further information on these topics are included in the root cause analysis for reference.

#### *5.3.3 Results of the collaborators retained for interpretation*

The results obtained with 10 collaborators kept for interpretation are presented in Table 20 (reference method) and Table 21 (alternative method).

**Table 20 - Positive results by the reference method (without participant 4 and 7)**

| Participant | Contamination level |     |     |
|-------------|---------------------|-----|-----|
|             | L0                  | L1  | L2  |
| 2           | 0/8                 | 6/8 | 8/8 |
| 3A          | 2/8                 | 6/8 | 8/8 |
| 3B          | 1/8                 | 7/8 | 8/8 |
| 5A          | 3/8                 | 6/8 | 8/8 |

| Participant  | Contamination level         |                              |                              |
|--------------|-----------------------------|------------------------------|------------------------------|
|              | L0                          | L1                           | L2                           |
| 5B           | 2/8                         | 6/8                          | 8/8                          |
| 6A           | 0/8                         | 6/8                          | 8/8                          |
| 6B           | 0/8                         | 6/8                          | 8/8                          |
| 8A           | 0/8                         | 6/8                          | 8/8                          |
| 8B           | 0/8                         | 6/8                          | 8/8                          |
| 10           | 0/8                         | 6/8                          | 8/8                          |
| <b>TOTAL</b> | <b>P<sub>0</sub> = 8/80</b> | <b>P<sub>1</sub> = 60/80</b> | <b>P<sub>2</sub> = 80/80</b> |

**Table 21 - Positive results (before and after confirmation) by the alternative methods (without participant 4 and 7)**

| Collaborator | Contamination level        |                           |                            |                           |                            |                           |
|--------------|----------------------------|---------------------------|----------------------------|---------------------------|----------------------------|---------------------------|
|              | L0                         |                           | L1                         |                           | L2                         |                           |
|              | <i>Before confirmation</i> | <i>After confirmation</i> | <i>Before confirmation</i> | <i>After confirmation</i> | <i>Before confirmation</i> | <i>After confirmation</i> |
| 2            | 0/8                        | 0/8                       | 6/8                        | 6/8                       | 8/8                        | 8/8                       |
| 3A           | 0/8                        | 0/8                       | 6/8                        | 6/8                       | 8/8                        | 8/8                       |
| 3B           | 0/8                        | 0/8                       | 6/8                        | 6/8                       | 8/8                        | 8/8                       |
| 5A           | 0/8                        | 0/8                       | 6/8                        | 6/8                       | 8/8                        | 8/8                       |
| 5B           | 0/8                        | 0/8                       | 6/8                        | 6/8                       | 8/8                        | 8/8                       |
| 6A           | 0/8                        | 0/8                       | 6/8                        | 6/8                       | 8/8                        | 8/8                       |
| 6B           | 0/8                        | 0/8                       | 6/8                        | 6/8                       | 8/8                        | 8/8                       |
| 8A           | 0/8                        | 0/8                       | 6/8                        | 6/8                       | 8/8                        | 8/8                       |

| Collaborator | Contamination level        |                             |                              |                               |                              |                               |
|--------------|----------------------------|-----------------------------|------------------------------|-------------------------------|------------------------------|-------------------------------|
|              | L0                         |                             | L1                           |                               | L2                           |                               |
|              | Before confirmation        | After confirmation          | Before confirmation          | After confirmation            | Before confirmation          | After confirmation            |
| 8B           | 0/8                        | 0/8                         | 6/8                          | 6/8                           | 8/8                          | 8/8                           |
| 10           | 0/8                        | 0/8                         | 6/8                          | 6/8                           | 8/8                          | 8/8                           |
| <b>TOTAL</b> | <b>P<sub>0</sub> = 0/8</b> | <b>CP<sub>0</sub> = 0/8</b> | <b>P<sub>1</sub> = 60/80</b> | <b>CP<sub>1</sub> = 60/80</b> | <b>P<sub>2</sub> = 80/80</b> | <b>CP<sub>2</sub> = 80/80</b> |

#### 5.3.4 Calculation of the specificity percentage (SP)

The percentage specificities (SP) of the reference method and of the alternative method, using the data after confirmation, based on the results of level L0 are the following (See Table 22).

**Table 22 - Percentage specificity**

|   |   |      |
|---|---|------|
| <b>Specificity for the reference method</b>   | $SP_{ref} = \left( 1 - \left( \frac{P_0}{N_-} \right) \right) \times 100 \% =$  | 90%  |
| <b>Specificity for the alternative method</b> | $SP_{alt} = \left( 1 - \left( \frac{CP_0}{N_-} \right) \right) \times 100 \% =$ | 100% |

N - number of all L0 tests

P<sub>0</sub> - total number of false-positive results obtained with the blank samples before confirmation

CP<sub>0</sub> - total number of false-positive results obtained with the blank samples

#### 5.3.5 Calculation of the sensitivity ( $SE_{alt}$ ), the sensitivity for the reference method ( $SE_{ref}$ ), the relative trueness (RT) and the false positive ratio for the alternative method (FPR)

Fractional positive results were obtained for the low inoculation level (L1).

A summary of the results of the collaborators retained for interpretation and obtained with the reference and the alternative methods for Level 1 is provided in Table 23.



**Table 23 - Summary of the obtained results with the reference method and the alternative method for Level 1 and/or Level 2**

| Level | Response                         | Reference method positive (R+)            | Reference method negative (R-)            |
|-------|----------------------------------|---|---|
| 1     | Alternative method positive (A+) | Positive agreement (A+/R+) <b>PA = 60</b> | Positive deviation (R-/A+) <b>PD = 0</b>  |
|       | Alternative method negative (A-) | Negative deviation (A-/R+) <b>ND = 1</b>  | Negative agreement (A-/R-) <b>NA = 19</b> |

Based on the data summarized in Table 10, the values of sensitivity of the alternative and reference methods, as well as the relative trueness and false positive ratio for the alternative method taking account the confirmations, are the following (See Table 24).

**Table 24 - Sensitivity, relative trueness and false positive ratio percentages**

|  |  | Level 1 |
|--|--|---------|
| <b>Sensitivity for the alternative method:</b>         | $SE_{alt} = \frac{(PA+PD)}{(PA+PD+ND)} \times 100\% =$ | 98.4    |
| <b>Sensitivity for the reference method:</b>           | $SE_{ref} = \frac{(PA+ND)}{(PA+PD+ND)} \times 100\% =$ | 100     |
| <b>Relative trueness</b>                               | $RT = \frac{(PA+NA)}{N} \times 100\% =$                | 98.75   |
| <b>False positive ratio for the alternative method</b> | $FPR = \frac{FP}{NA} \times 100\% =$                   | 0       |

### 5.3.6 Interpretation of data

The negative deviations are listed in Table 25.

**Table 25 - Negative deviations for Level 1**

| Sample n°      | Reference method result | Alternative method results | Inoculation (CFU/ml) |
|----------------|-------------------------|----------------------------|----------------------|
| Participant 3B | +                       | ND                         | 11                   |

For an **unpaired study design**, the difference between (ND – PD) is calculated for the level(s) where fractional recovery is obtained (so  $L_1$  and possibly  $L_2$ ). The observed value found for (ND – PD) shall not be higher than the AL. The AL is defined as  $[(ND - PD)_{max}]$  and calculated per level where fractional recovery is obtained as described below using the following three parameters:

$$(p+)_{ref} = \frac{P_x}{N_x}$$

where

$P_x$  = number of samples with a positive result obtained with the reference method at level  $x$  ( $L_1$  or  $L_2$ ) for all the collaborators

$N_x$  = number of samples tested at level  $x$  ( $L_1$  or  $L_2$ ) with the reference method by all the collaborators

$$(p+)_{alt} = \frac{CP_x}{N_x}$$

where

$CP_x$  = number of samples with a confirmed positive result obtained with the alternative method at level  $x$  ( $L_1$  or  $L_2$ ) for all the collaborators;

$N_x$  = number of samples tested at level  $x$  ( $L_1$  or  $L_2$ ) with the alternative method by all the collaborators.

$$(ND-PD)_{max} = \sqrt{3N_x \times \left( (p+)_{ref} + (p+)_{alt} - 2((p+)_{ref} \times (p+)_{alt}) \right)}$$

where

$N_x$  = number of samples tested for level  $x$  ( $L_1$  or  $L_2$ ) with the reference method by all the collaborators.

In this study, fractional recovery was observed at Level 1. The calculations, according to the EN ISO 16140-2:2016, are shown in Table 26.

**Table 26 - Calculations**

|                       | Level 1     |
|-----------------------|-------------|
| $N_x$                 | 80          |
| $(p+)_{ref}$          | 0.76        |
| $(p+)_{alt}$          | 0.75        |
| $AL = (ND - PD) \max$ | 18          |
| $ND - PD$             | 1           |
| <b>Conclusion</b>     | <b>Pass</b> |

The ISO 16140-2 (2016) requirements are fulfilled as  $(ND - PD)$  is below the Acceptability Limit.

#### 5.3.7 Evaluation of the RLOD between laboratories

The RLOD was calculated using the EN ISO 16140-2:2016 Excel spreadsheet available at <http://standards.iso.org/iso/16140> - RLOD (clause 5-1-4-2 Calculation and interpretation of RLOD) version 06.07.2015. The results are used only for information (see Table 1).

**Table 1 - RLOD**

| Results by method  |   |                            |                                     |  |                   |                   |  |                   |                   |
|--------------------|---|----------------------------|-------------------------------------|--|-------------------|-------------------|--|-------------------|-------------------|
| Method             | Method effect<br>$F_i$  | Log method effect<br>$f_i$ | SD of log method effect<br>$s_{fi}$ | LOD <sub>50%</sub> = 50% limit of detection in cfu per sample size |                   |                   | LOD <sub>95%</sub> = 95% limit of detection in cfu per sample size |                   |                   |
|                    |   |                            |                                     | Detection limit  | Lower conf. limit | Upper conf. limit | Detection limit  | Lower conf. limit | Upper conf. limit |
|                    |   |                            |                                     | $d_{0.5,i}$  | $d_{0.5,i,L}$     | $d_{0.5,i,U}$     | $d_{0.95,i}$   | $d_{0.95,i,L}$    | $d_{0.95,i,U}$    |
| Reference          | 0.143   | -1.945                     | 0.118                               | 4.85   | 3.83              | 6.13              | 20.95  | 16.56             | 26.51             |
| Alternative        | 0.139   | -1.973                     | 0.117                               | 4.98   | 3.95              | 6.29              | 21.54  | 17.06             | 27.20             |
| <b>Conclusions</b> | The methods are <u>not significantly</u> different at the 5% significance level (change in deviance of the model with method effects to the null model $D_{method} = 0.03$ with 1 degree of freedom, p-value 0.86). |                            |                                     |  |                   |                   |  |                   |                   |
|                    | The relative limit of detection (RLOD) of the alternative method, as compared to the reference method, is 1.03 with a 90% confidence interval of 0.78 - 1.35.   |                            |                                     |  |                   |                   |  |                   |                   |

#### 5.4 Root cause analysis – interlaboratory study

A root cause analysis was carried out to investigate issues that occurred during the ILS. The two issues included in the root cause analysis were the number of external study participants and the positive blank samples obtained with the reference method.

#### 5.4.1 Recruitment of Participants

There was difficulty in recruiting the required number of participants for the ILS due to the nature of the method (i.e. instrument-based). The study was arranged on two previous occasions however they were cancelled due to a lack of participants.

On this occasion, 14 external participants were recruited for the study included in the current report. The parcels for 13 out of 14 participants arrived in time and 2 participants did not complete the testing. As a result, 11 external datasets were received from the 14 participants recruited for the ILS.

From the 11 external datasets that were received, there were 2 datasets that required further investigation.

#### Participant 7

No detection time was recorded with the Soleris instrument for any samples analysed in the ILS. However, the vials were streaked, and fractional positive confirmation results were received for level 1. The vials also showed visible signs of growth (e.g. change of plug colour, turbidity in broth, as shown in Figure 1).

Further investigation found that the graphs were not plotting, indicating that the instrument was not recording measurements. The instrument was out of calibration and was using a previous software version. All other participants in the study were using the same software version.

Due to the instrument issues, the results for participant 7 were excluded from the analysis.




Figure 1 - Image taken by participant showing (left) positive vial and (right) negative vial

#### Participant 4

It was noted that participant 4 stored their samples frozen for 2 weeks instead of 1 week stated in the ILS instructions. The impact of the increased storage time was further investigated below

The plan for preparation, storage and testing of the samples was as follows:

| Monday 10 <sup>th</sup> June         | Tuesday 11 <sup>th</sup> June     | Wednesday 12 <sup>th</sup> June –<br>Thursday 13 <sup>th</sup> June                           | Monday 17 <sup>th</sup> June  |
|--------------------------------------|-----------------------------------|---|---|
| Samples were inoculated, then frozen | Samples were packaged and shipped | Samples arrived in participating laboratories. Participants froze the samples until Analysis. | Participants defrosted under controlled conditions and began testing. |

Participant 7 stored the samples for an additional week before testing on Monday 24<sup>th</sup> June which was beyond the scope of the stability trial. The storage had a potential impact on reference method, the number of positives for each method is shown below:

- L1 2/8 ref, 6/8 alt
- L2 4/8 ref, 8/8 alt

In addition to the potential stability issue, there were 5 positive blank samples for the reference method, and 1 positive blank sample for the alternative method. The participant provided images of the reference plates and alternative confirmation plates. The images were analysed by the expert laboratory for colony morphology. Figure 2 shows the positive blank sample for the alternative method, which has similar morphology to the inoculating strain (*Staphylococcus saprophyticus*). The investigation carried out on the images provided by the participant were inconclusive, however it indicates that there could have been contamination between inoculated and uninoculated samples.



*Figure 2 – Alternative confirmation plate for Participant 7, sample 10 which shows contamination of a blank sample with colony morphology similar to that of the inoculating strain*

Due to the results of the investigations and the potential impact of additional storage on the reference method, the results for participant 4 have been excluded.

#### Use of an internal participant

Due to the reduction in datasets received compared to the number of participants that were recruited and the exclusion of 2 datasets, there were 9 external datasets remaining for analysis. The minimum number of datasets required for the analysis is 10.

In addition to 14 external participants, an internal participant from the analytical services team at Campden BRI, also took part in the study. This participant was in addition to the expert laboratory participant. The internal participant was not involved in the preparation of samples and had not performed the method before the interlaboratory study therefore could be considered to be an independent participant.

The decision-making process to include the internal dataset for the analysis was guided by the relevant definitions within ISO16140-1 (2016) and ISO 16140-2 (2016) concerning the organisation of the ILS and the participants taking part in the study.

In section 2.13 of ISO16140-1 (2016), a collaborator is defined as an individual laboratory technician who works completely independently for the other collaborators using a different set of blind samples or test portions.

The organising laboratory is defined in section 2.45 of ISO16140-1 (2106), as an expert independent laboratory with the responsibility for managing all of the technical and statistical analysis involved in the validation study including the method comparison study and interlaboratory study.

Section 5.2.2 of ISO 16140-2 (2016) outlines the measurement protocol of the ILS, which states that that the technicians involved in the preparation of the samples used in the ILS shall not take part in the testing of those samples within the interlaboratory study

During the study, the samples analysed at the organising laboratory were tested by a technician who had not been involved in the preparation of the samples. As this set of samples were analysed by a technician who was not responsible for the preparation of the samples or the data or statistical analysis of the samples, it is proposed that this data set could be considered as independent from the other collaborators.

The statistical analysis was performed on 9 external datasets and 1 internal dataset performed by analytical services.

The expert laboratory dataset was not included in the analysis.

#### *5.4.2 Positive blank samples*




Of the 10 datasets included in the analysis, 4 participants had up to 3 blank samples with a positive result. 8 out of 80 uninoculated samples were positive for the reference method only. The number of positive blank samples per participant is shown in Table 28.

**Table 28 – Number of positive uninoculated samples obtained by each participant**




| Participant  | Contamination level         |
|--------------|-----------------------------|
|              | L0                          |
| 2            | 0/8                         |
| 3A           | 2/8                         |
| 3B           | 1/8                         |
| 5A           | 3/8                         |
| 5B           | 2/8                         |
| 6A           | 0/8                         |
| 6B           | 0/8                         |
| 8A           | 0/8                         |
| 8B           | 0/8                         |
| 10           | 0/8                         |
| <b>TOTAL</b> | <b>P<sub>0</sub> = 8/80</b> |



Further investigations were carried out for participants with positive blank samples. Participants with positive blank samples were asked to take images of plates and carry out MALDI-ToF confirmations, if available. The results of the investigations are shown in Table 29.

**Table 29 – Investigation of positive blank samples**

| Participant | Sample | Count on plate | Image of reference plate taken by the participant                                    | MALDI-ToF result  |
|-------------|--------|----------------|--|---|
| 3A          | 13     | 2              |    | Colony A - <i>Bacillus cereus</i> 2.23<br>Colony B – <i>Corynebacterium mucifaciens</i> 1.98 (lower confidence score) |
|             | 22     | 2              |   | Colony A – <i>Bacillus cereus</i> 2.3<br>Colony B – <i>Bacillus cereus</i> 2.3  |
| 3B          | 1      | 2              |  | Colony A – <i>Bacillus cereus</i> 2.3<br>Colony B – <i>Bacillus cereus</i> 2.11                                       |



|    |    |   |  |     |
|----|----|---|--|-----|
| 5A | 1  | 1 |    | N/A |
|    | 13 | 1 |   | N/A |
|    | 16 | 1 |  | N/A |

|    |    |               |   |     |
|----|----|---------------|---|-----|
| 5B | 7  | 11<br>(small) |   | N/A |
|    | 10 | 2             |  | N/A |

The investigations showed minimal contamination on each plate for the reference method. The colony morphologies/identifications obtained indicate that these samples were not contaminated with the inoculating strain. In conclusion, the blank samples with non target organisms present were retained for statistical analysis as they were considered to be a not detected result.

#### Root cause analysis – positive blank samples

The scope of the alternative method is the detection of total viable count. Positive blank samples were obtained with the reference method only. For each sample, 1ml of the 1 in 10 dilution was plated on 1x petri dish for the reference method.

Some potential causes for the positive blank samples are listed in Table 30.

**Table 30 – Potential root causes for positive blank samples**

| Potential Cause                            | Reasons for/against potential cause  | Likelihood of potential cause to produce positive blank samples |
|--|--|---|
| Competency/familiarity of reference method | <p>Participants were asked to complete a survey on how frequently they perform the reference method.</p> <p>2/4 affected participants have completed the survey and responded that they are performing the reference method on a daily/weekly basis.</p> <p>The results of the remaining participants will be presented at the MVTC.</p> | Unlikely  |
| Sterility of media                         | <p>1 participant found 1 bottle of PCA that was visibly contaminated. The contaminated bottle was not used in the study, although there is possibility that more bottles could have been contaminated.</p> <p>The occurrence of positive blank samples was limited to 2 laboratories.</p>  | Likely  |
| Handling of samples                        | <p>Images/confirmations show that the contamination observed in the positive blank samples is different to the inoculating strain. This indicates that inoculated and uninoculated samples have not been contaminated with each other.</p>   | Unlikely  |
| Aseptic technique                          | <p>Images show colonies on the surface of the plate, this indicates that the samples could have been contaminated after the sample was plated e.g. when pouring media.</p> <p>The occurrence of positive blank samples was limited to 2 laboratories.</p>  | Possible – when pouring media                                   |

|                                   |   |          |
|-----------------------------------|---|----------|
| Sample was naturally contaminated | <p>The sample selected is a UHT product, and therefore less likely to be naturally contaminated.</p> <p>Natural contamination was not observed in screening and the stability study.</p> <p>Samples were prepared to minimise potential contamination:</p> <ul style="list-style-type: none"> <li>• Weighed in a laminar flow</li> <li>• Uninoculated samples were frozen immediately after weighing</li> <li>• Inoculated samples were weighed on the day of inoculation</li> </ul> <p>Images from the participants indicate that contamination was typically on the surface of the plate, this does not indicate natural contamination of the sample.</p> | Unlikely |
|-----------------------------------|---|----------|

Two other validation studies involving the detection of low levels of total viable count have been completed, 2021LR100 NF105 Commercial Sterility vials and 2021LR94 Certablue TVC vials. False positive blank samples were also found during the interlaboratory study for both validation studies. Both of these studies were accepted by the MVTC following review of the results and potential causes of the positive blank samples. For the Certablue TVC vials, the contamination was found in the alternative method only. As a result of discussions at the MVTC, a warning was added to the kit insert to emphasize the risk of cross contamination.

#### 5.4.3 Summary

14 external participants were recruited for the current ILS and 11 external datasets were received. After initial review, 2 external datasets were excluded and the statistical analysis was carried out with 9 external dataset + 1 internal dataset. The internal dataset was carried out by a member of the analytical services team who was not involved in the validation study and had not previously used the NF TVC vial method. In this case it was considered that the internal dataset included in the study was independent.

Positive samples were obtained by the reference method in 10% of blank samples tested. The contaminated samples were limited to 4 participants in 2 laboratories. Further investigations and a root cause analysis were carried out, which identified that the contamination was likely due to the sterility of PCA used or aseptic technique. The impact of this on the validity of the results is minimal.

In conclusion, the interlaboratory study data and interpretations comply with the ISO 16140-2:2016 requirements.

## 6 CONCLUSION

The method comparison study conclusions are:

- The observed values for ND-PD for the individual categories and for all categories meet the acceptability limits (observed values  $\leq$  AL).
- The RLOD values (using the confirmed alternative method results) meet the acceptability limit, which is 2.5 for unpaired studies, for all categories tested.
- The alternative Neogen Soleris® NF-TVC detection method is selective and specific.

The **inter-laboratory study conclusions** are:

- The observed value for ND-PD is lower than the acceptability limits.
- There are no individual categories tested for the IL.
- The data and interpretations comply with the EN ISO 16140-2:2016 requirements.

**The Neogen NF-TVC vial is considered equivalent to the ISO standard.**


Date,

17/10/2024


Signature,

Alice Foxall – Project Manager, Molecular Microbiology and Methods

**ANNEX A: Flow diagram of the reference method**



## ANNEX B: Flow diagram of the alternative method



Software settings

| Threshold | Skip | Shuteye | Temperature | Test duration |
|-----------|------|---------|-------------|---------------|
| 10        | 1    | 30      | 30°C        | 30 hours      |

The confirmation step is **not required** for the user, and was performed for purposes of this certification only



### ANNEX C: Raw data from the sensitivity study

| Milk and dairy products |                  |   |                   |                |   |                  |           |
|-------------------------|------------------|---|-------------------|----------------|---|------------------|-----------|
| Sample code             | Type             | Sample used                                   | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| T61                     | Pasteurised milk | Semi skimmed milk 1                           | -6                | -              | -   | -                | NA        |
| T62                     | Pasteurised milk | Semi skimmed milk 2                           | -6                | -              | -   | -                | NA        |
| T63                     | Pasteurised milk | Jersey milk                                   | -6                | +              | +   | +                | PA        |
| T64                     | Pasteurised milk | Whole milk                                    | -6                | -              | -   | -                | NA        |
| T65                     | Pasteurised milk | Skimmed milk                                  | -6                | +              | +   | +                | PA        |
| D6                      | Pasteurised milk | Semi skimmed milk that tastes like whole      | -6                | -              | -   | -                | NA        |
| D7                      | Pasteurised milk | West country whole milk                       | -6                | +              | +   | +                | PA        |
| D8                      | Pasteurised milk | Organic skimmed milk                          | -6                | -              | -   | -                | NA        |
| D9                      | Pasteurised milk | Jersey milk taste the difference              | -6                | +              | +   | +                | PA        |
| D10                     | Pasteurised milk | Organic fresh semi skimmed milk               | -6                | -              | -   | -                | NA        |
| D11                     | Pasteurised milk | Fresh whole milk fresher for longer           | -6                | -              | -   | -                | NA        |
| D12                     | Pasteurised milk | Fresh skimmed milk fresher for longer         | -6                | -              | -   | -                | NA        |
| D13                     | Pasteurised milk | British whole milk so organic                 | -6                | +              | +   | +                | PA        |
| D14                     | Pasteurised milk | Skimmed milk tastes like semi skimmed         | -6                | -              | -   | -                | NA        |
| D15                     | Pasteurised milk | British semi skimmed milk so organic          | -6                | +              | +   | +                | PA        |
| D16                     | Pasteurised milk | Fresh semi skimmed milk 2l fresher for longer | -6                | -              | -   | -                | NA        |
| D17                     | Pasteurised milk | British skimmed milk                          | -6                | -              | -   | -                | NA        |
| D18                     | Pasteurised milk | British whole milk                            | -6                | -              | -   | -                | NA        |
| D19                     | Pasteurised milk | British semi skimmed milk                     | -6                | -              | -   | -                | NA        |
| D20                     | Pasteurised milk | British whole milk                            | -6                | +              | +   | +                | PA        |





| Milk and dairy products |                            |  |                   |                |   |                  |           |
|-------------------------|----------------------------|--|-------------------|----------------|---|------------------|-----------|
| Sample code             | Type                       | Sample used                              | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| T66                     | Pasteurised dairy products | Fresh whipping cream                     | -4                | -              | -   | -                | NA        |
| T67                     | Pasteurised dairy products | Soured cream                             | -4                | +              | +   | -                | PD        |
| T68                     | Pasteurised dairy products | Fat free cottage cheese, onion & chive   | -4                | +              | +   | +                | PA        |
| T69                     | Pasteurised dairy products | Butter milk                              | -4                | +              | +   | +                | PA        |
| T70                     | Pasteurised dairy products | Light original soft cheese               | -4                | -              | -   | -                | NA        |
| D26                     | Pasteurised dairy products | Chocolate milkshake                      | -4                | +              | +   | +                | PA        |
| D27                     | Pasteurised dairy products | Strawberry milkshake                     | -4                | -              | -   | -                | NA        |
| D28                     | Pasteurised dairy products | Milk chocolate chip ice crea             | -4                | +              | +   | +                | PA        |
| D29                     | Pasteurised dairy products | Caramel and vanilla ice cream            | -4                | -              | -   | -                | NA        |
| D30                     | Pasteurised dairy products | Vanilla ice cream                        | -4                | -              | -   | -                | NA        |
| D31                     | Pasteurised dairy products | Ricotta cheese                           | -4                | -              | -   | -                | NA        |
| D32                     | Pasteurised dairy products | Mozzarella cheese                        | -4                | +              | +   | +                | PA        |
| D33                     | Pasteurised dairy products | Garlic and herb soft french cream cheese | -4                | -              | -   | -                | NA        |
| D34                     | Pasteurised dairy products | Liquid cheese original                   | -4                | -              | -   | -                | NA        |
| D35                     | Pasteurised dairy products | Organic double free range cream          | -4                | +              | +   | +                | PA        |
| D36                     | Pasteurised dairy products | Jersey double cream taste the difference | -4                | +              | +   | +                | PA        |
| D37                     | Pasteurised dairy products | Double organic cream                     | -4                | +              | +   | +                | PA        |
| D38                     | Pasteurised dairy products | Cream                                    | -4                | -              | -   | -                | NA        |
| D39                     | Pasteurised dairy products | British double cream                     | -4                | -              | -   | -                | NA        |
| D40                     | Pasteurised dairy products | Clotted cream                            | -4                | -              | -   | -                | NA        |
| T71                     | Dry milk products          | Dried skim milk                          | -1                | -              | -   | -                | NA        |
| T72                     | Dry milk products          | Milk shake powder                        | -1                | +              | +   | +                | PA        |



| Milk and dairy products |                   |                                       |                   |                |   |                  |           |
|-------------------------|-------------------|---------------------------------------|-------------------|----------------|---|------------------|-----------|
| Sample code             | Type              | Sample used                           | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| T73                     | Dry milk products | Skimmed milk powder                   | -1                | -              | -   | -                | NA        |
| T74                     | Dry milk products | Non-fat skim milk                     | -1                | +              | +   | +                | PA        |
| T75                     | Dry milk products | Non-fat skim milk                     | -1                | +              | +   | +                | PA        |
| D46                     | Dry milk products | Strawberry flavoured milk powder      | -1                | -              | -   | -                | NA        |
| D47                     | Dry milk products | Vanilla flavoured milk powder         | -1                | -              | -   | -                | NA        |
| D48                     | Dry milk products | Banana milkshake powder tub           | -1                | -              | -   | -                | NA        |
| D49                     | Dry milk products | Chocolate milkshake powder tub        | -1                | +              | +   | +                | PA        |
| D50                     | Dry milk products | Strawberry milkshake powder tub       | -1                | -              | -   | -                | NA        |
| D51                     | Dry milk products | Chocolate flavour dessert             | -1                | +              | +   | +                | PA        |
| D52                     | Dry milk products | Butterscotch flavour dessert          | -1                | +              | +   | +                | PA        |
| D53                     | Dry milk products | Strawberry flavour dessert            | -1                | +              | +   | +                | PA        |
| D54                     | Dry milk products | Custard powder                        | -1                | +              | +   | +                | PA        |
| D55                     | Dry milk products | Instant custard powder no added sugar | -1                | +              | +   | +                | PA        |
| D56                     | Dry milk products | Instant custard powder 2              | -1                | +              | +   | +                | PA        |
| D57                     | Dry milk products | Instant custard powder 3              | -1                | +              | +   | +                | PA        |
| D58                     | Dry milk products | Strawberry milk powder                | -1                | +              | +   | +                | PA        |
| D59                     | Dry milk products | Chocolate milk powder                 | -1                | +              | +   | +                | PA        |
| D60                     | Dry milk products | Speculoos whey powder                 | -1                | +              | +   | +                | PA        |



| Fishery products  |                 |  |                   |                |   |                  |           |
|---|-----------------|--|-------------------|----------------|---|------------------|-----------|
| Samples in bold required 72h incubation for confirmation result |                 |  |                   |                |   |                  |           |
| Sample code   | Type            | Sample used                              | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| S41   | Raw fish        | Salmon fillet                            | -8                | +              | +   | -                | PD        |
| S42   | Raw fish        | Scottish tail fillets                    | -8                | -              | -   | +                | ND        |
| S43   | Raw fish        | Yellowfin Tuna Steaks                    | -8                | -              | -   | -                | NA        |
| S44   | Raw fish        | Icelandic cod fillets                    | -8                | -              | -   | -                | NA        |
| S45   | Raw fish        | Scottish Loch trout fillets              | -8                | -              | -   | -                | NA        |
| T16   | Raw fish        | Basa fillets                             | -8                | +              | +   | +                | PA        |
| T17   | Raw fish        | Sea bass fillets                         | -8                | -              | -   | +                | ND        |
| <b>T18</b>  | <b>Raw fish</b> | <b>Smoked basa fillets</b>               | <b>-8</b>         | <b>+</b>       | <b>+</b>                                  | <b>+</b>         | <b>PA</b> |
| T19   | Raw fish        | Smoked Norwegian haddock fillets         | -8                | +              | +   | +                | PA        |
| T20   | Raw fish        | Salmon fillets                           | -8                | +              | +   | +                | PA        |
| F61   | Raw fish        | Wild salmon fillets                      | -8                | -              | -   | -                | NA        |
| F62   | Raw fish        | Raw squid rings                          | -8                | -              | -   | +                | ND        |
| <b>F63</b>  | <b>Raw fish</b> | <b>Cod fillets</b>                       | <b>-8</b>         | <b>+</b>       | <b>+</b>                                  | <b>+</b>         | <b>PA</b> |
| F64   | Raw fish        | White fish fillets                       | -8                | +              | +   | +                | PA        |
| <b>F65</b>  | <b>Raw fish</b> | <b>Boneless Seabass Fillets</b>          | <b>-8</b>         | <b>+</b>       | <b>+</b>                                  | <b>+</b>         | <b>PA</b> |
| <b>F66</b>  | <b>Raw fish</b> | <b>Mussels</b>                           | <b>-8</b>         | <b>+</b>       | <b>+</b>                                  | <b>+</b>         | <b>PA</b> |
| F67   | Raw fish        | Salmon tails                             | -8                | +              | +   | +                | PA        |
| F68   | Raw fish        | Scottish salmon fillets                  | -8                | -              | -   | -                | NA        |
| F69   | Raw fish        | Mild and delicate Scottish smoked salmon | -8                | -              | -   | -                | NA        |



| Fishery products  |                 |   |                   |                |   |                  |           |
|---|-----------------|---|-------------------|----------------|---|------------------|-----------|
| Samples in bold required 72h incubation for confirmation result |                 |   |                   |                |   |                  |           |
| Sample code   | Type            | Sample used   | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| <b>F70</b>  | <b>Raw fish</b> | <b>Skin on cod fillets</b>                                  | <b>-8</b>         | +              | +   | +                | PA        |
| <b>F16</b>  | <b>Raw fish</b> | <b>Cod fillets 2</b>  | <b>-8</b>         | +              | +   | +                | PA        |
| F17   | Raw fish        | Smoked haddock fillets                                      | -8                | +              | +   | +                | PA        |
| F18   | Raw fish        | Wild salmon fillets   | -8                | -              | -   | -                | NA        |
| F19   | Raw fish        | Norwegian Flaky cod fillets                                 | -8                | +              | -   | -                | PPNA      |
| F20   | Raw fish        | White fish fillets  | -8                | -              | -   | -                | NA        |
| S61   | RTC fish        | Smoked kipper fillets w/ butter                             | -8                | -              | -   | -                | NA        |
| S62   | RTC fish        | Mediterranean tomato fish bakes                             | -8                | -              | -   | -                | NA        |
| S63   | RTC fish        | Cod & spinach gratin  | -8                | -              | -   | -                | NA        |
| <b>S64</b>  | <b>RTC fish</b> | <b>Smoked haddock fishcakes w/ cheddar &amp; leek sauce</b> | <b>-8</b>         | +              | +   | +                | PA        |
| S65   | RTC fish        | Chunky cod fish fingers                                     | -8                | -              | -   | -                | NA        |
| T21   | RTC fish        | King prawn flavour tails                                    | -8                | -              | -   | -                | NA        |
| T22   | RTC fish        | Breaded Cod Fishcakes                                       | -8                | -              | -   | -                | NA        |
| T23   | RTC fish        | Cod fishcakes with roasted tomato and mozzarella            | -8                | -              | -   | -                | NA        |
| <b>T24</b>  | <b>RTC fish</b> | <b>Golden breaded cod fillets</b>                           | <b>-8</b>         | +              | +   | +                | PA        |
| T25   | RTC fish        | Smoked haddock fishcakes with cheddar and leek              | -8                | -              | -   | -                | NA        |
| <b>F6</b>   | <b>RTC fish</b> | <b>Salmon fishcakes</b>                                     | <b>-8</b>         | +              | +   | +                | PA        |
| F7  | RTC fish        | Creamy cod fishcakes  | -8                | -              | -   | +                | ND        |
| F11   | RTC fish        | Crisp lightly dusted Basa fillets                           | -8                | +              | +   | +                | PA        |
| F13   | RTC fish        | Smoked haddock fishcakes with melting cheddar centre        | -8                | -              | -   | +                | ND        |
| F14   | RTC fish        | Breaded white fish fillets                                  | -8                | -              | -   | -                | NA        |
| F15   | RTC fish        | Salmon fishcakes with melting hollandaise centre            | -8                | +              | +   | +                | PA        |
| <b>F32</b>  | <b>RTC fish</b> | <b>Cod &amp; Chorizo Fishcakes</b>                          | <b>-8</b>         | +              | +   | +                | PA        |
| F33   | RTC fish        | Hot Smoked Salmon Fillets                                   | -8                | -              | -   | -                | NA        |
| <b>F34</b>  | <b>RTC fish</b> | <b>Smoked Haddock Fishcakes</b>                             | <b>-8</b>         | +              | +   | +                | PA        |



| Fishery products  |                    |   |                   |                |   |                  |           |
|---|--------------------|---|-------------------|----------------|---|------------------|-----------|
| Samples in bold required 72h incubation for confirmation result |                    |   |                   |                |   |                  |           |
| Sample code   | Type               | Sample used                             | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| F35   | RTC fish           | Finest Prawn Cocktail                   | -8                | -              | -   | -                | NA        |
| F36   | RTC fish           | Salmon Spinach & Lemon Fishcakes        | -8                | -              | -   | -                | NA        |
| F37   | RTC fish           | Cod Fishcakes                           | -8                | +              | +   | +                | PA        |
| F38   | RTC fish           | Flipper Dippers                         | -8                | +              | +   | +                | PA        |
| F39   | RTC fish           | Honey Roast Salmon Flakes               | -8                | -              | -   | -                | NA        |
| F40   | RTC fish           | Thai Inspired Salmon & Cod Fishcakes    | -8                | +              | +   | +                | PA        |
| <b>T26</b>  | <b>Crustaceans</b> | <b>Cold water large prawns</b>          | <b>-6</b>         | <b>+</b>       | <b>+</b>                                  | <b>-</b>         | <b>PD</b> |
| S1  | Crustaceans        | White Crab                              | -6                | +              | +   | -                | PD        |
| T28   | Crustaceans        | Garlic & herb king prawns               | -6                | -              | -   | -                | NA        |
| T29   | Crustaceans        | Soy, ginger & chilli king prawns        | -6                | -              | -   | -                | NA        |
| T30   | Crustaceans        | Raw & peeled jumbo king prawns          | -6                | -              | -   | -                | NA        |
| S2  | Crustaceans        | Brown Crab                              | -6                | -              | -   | -                | NA        |
| F10   | Crustaceans        | Cold water large prawns                 | -6                | -              | -   | -                | NA        |
| S3  | Crustaceans        | Crab                                    | -6                | -              | -   | -                | NA        |
| S4  | Crustaceans        | Premium Crab meat                       | -6                | +              | +   | -                | PD        |
| S5  | Crustaceans        | White Crab 2                            | -6                | -              | -   | +                | ND        |
| S6  | Crustaceans        | White & Brown Crab                      | -6                | -              | -   | +                | ND        |
| S7  | Crustaceans        | Crayfish tails                          | -6                | -              | -   | -                | NA        |
| S8  | Crustaceans        | Crab meat shredded                      | -6                | -              | -   | -                | NA        |
| F54   | Crustaceans        | Raw Chilli & Lime Marinated King Prawns | -6                | +              | +   | +                | PA        |
| F55   | Crustaceans        | Raw Headless King Prawns                | -6                | -              | -   | -                | NA        |
| S9  | Crustaceans        | Canadian cooked lobster                 | -6                | -              | -   | -                | NA        |
| F57   | Crustaceans        | Lemon & Pepper King Prawns              | -6                | +              | +   | +                | PA        |
| S10   | Crustaceans        | Organic Raw King Prawns                 | -6                | -              | -   | -                | NA        |
| F59   | Crustaceans        | Raw King Prawns                         | -6                | +              | +   | +                | PA        |

2022LR113 Neogen Soleris Non-Fermenting Total  
Viable Count (NF-TVC) vial MCS ILS Summary  
Report 14/10/2024



| Fishery products  |             |                         |                   |                |   |                  |           |
|---|-------------|-------------------------|-------------------|----------------|---|------------------|-----------|
| Samples in bold required 72h incubation for confirmation result |             |                         |                   |                |   |                  |           |
| Sample code   | Type        | Sample used             | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| F60   | Crustaceans | Raw & peeled red shrimp | -6                | +              | +   | +                | PA        |



| Fresh produce and fruits |  |                                      |                   |                |   |                  |           |
|--------------------------|--|--------------------------------------|-------------------|----------------|---|------------------|-----------|
| Sample code              | Type                                     | Sample used                          | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| T31                      | Cut ready-to-eat vegetables/leafy greens | Italian wild rocket                  | -8                | -              | -   | -                | NA        |
| T32                      | Cut ready-to-eat vegetables/leafy greens | Babyleaf salad                       | -8                | -              | -   | -                | NA        |
| T33                      | Cut ready-to-eat vegetables/leafy greens | Fine beans & tenderstem broccoli     | -8                | -              | -   | -                | NA        |
| T34                      | Cut ready-to-eat vegetables/leafy greens | Mixed leaf salad                     | -8                | -              | -   | -                | NA        |
| T35                      | Cut ready-to-eat vegetables/leafy greens | Butterhead salad                     | -8                | -              | -   | -                | NA        |
| P6                       | Cut ready-to-eat vegetables/leafy greens | Chopped spinach                      | -8                | +              | +   | +                | PA        |
| P7                       | Cut ready-to-eat vegetables/leafy greens | Edamame stir dry                     | -8                | +              | +   | +                | PA        |
| P8                       | Cut ready-to-eat vegetables/leafy greens | Baby vegetable stir dry              | -8                | +              | +   | +                | PA        |
| P9                       | Cut ready-to-eat vegetables/leafy greens | Fine beans and tenderstem broccoli   | -8                | +              | +   | +                | PA        |
| P10                      | Cut ready-to-eat vegetables/leafy greens | Mixed vegetable stir fry             | -8                | -              | -   | -                | NA        |
| P11                      | Cut ready-to-eat vegetables/leafy greens | Hot and spicy stir fry               | -8                | +              | +   | +                | PA        |
| P12                      | Cut ready-to-eat vegetables/leafy greens | Florette baby leaf & rocket          | -8                | +              | +   | +                | PA        |
| P13                      | Cut ready-to-eat vegetables/leafy greens | Mixed pepper stir fry veg mix        | -8                | -              | -   | -                | NA        |
| P14                      | Cut ready-to-eat vegetables/leafy greens | Tenderstem stir fry veg mix          | -8                | +              | +   | +                | PA        |
| P15                      | Cut ready-to-eat vegetables/leafy greens | Organic rocket salad                 | -8                | +              | +   | +                | PA        |
| P16                      | Cut ready-to-eat vegetables/leafy greens | Oriental vegetable stir dry          | -8                | +              | +   | +                | PA        |
| P17                      | Cut ready-to-eat vegetables/leafy greens | Wild rocket                          | -8                | +              | +   | +                | PA        |
| P18                      | Cut ready-to-eat vegetables/leafy greens | Spinach, watercress and rocket salad | -8                | +              | +   | +                | PA        |
| P19                      | Cut ready-to-eat vegetables/leafy greens | Cut babygem lettuce                  | -8                | +              | +   | +                | PA        |



| Fresh produce and fruits |  |  |                   |                |   |                  |           |
|--------------------------|--|--|-------------------|----------------|---|------------------|-----------|
| Sample code              | Type                                     | Sample used                            | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| P20                      | Cut ready-to-eat vegetables/leafy greens | Cut peppers and lettuce mix            | -8                | +              | +   | -                | PD        |
| T36                      | Fresh fruit/cut RTE fruit products       | Blueberries                            | -5                | -              | -   | -                | NA        |
| T37                      | Fresh fruit/cut RTE fruit products       | Strawberries                           | -5                | -              | -   | -                | NA        |
| T38                      | Fresh fruit/cut RTE fruit products       | Berry fruit salad                      | -5                | -              | -   | -                | NA        |
| T39                      | Fresh fruit/cut RTE fruit products       | Apple banana strawberry and grape      | -5                | -              | -   | -                | NA        |
| T40                      | Fresh fruit/cut RTE fruit products       | Melon kiwi and strawberry              | -5                | -              | -   | -                | NA        |
| P26                      | Fresh fruit/cut RTE fruit products       | Apples, mango strawberry and raspberry | -5                | +              | +   | +                | PA        |
| P27                      | Fresh fruit/cut RTE fruit products       | Mango                                  | -5                | +              | +   | +                | PA        |
| P28                      | Fresh fruit/cut RTE fruit products       | Pomegranate seeds                      | -5                | +              | +   | +                | PA        |
| P29                      | Fresh fruit/cut RTE fruit products       | Mango                                  | -5                | +              | +   | -                | PD        |
| P30                      | Fresh fruit/cut RTE fruit products       | Grapes and berries                     | -5                | +              | +   | +                | PA        |
| P31                      | Fresh fruit/cut RTE fruit products       | Rainbow fruit platter                  | -5                | +              | +   | +                | PA        |
| P32                      | Fresh fruit/cut RTE fruit products       | Watermelon fingers                     | -5                | +              | +   | +                | PA        |
| P33                      | Fresh fruit/cut RTE fruit products       | Melon stone fruit and raspberry        | -5                | +              | +   | +                | PA        |
| P34                      | Fresh fruit/cut RTE fruit products       | Fruit platter                          | -5                | +              | +   | +                | PA        |
| P35                      | Fresh fruit/cut RTE fruit products       | Melon and pineapple fingers            | -5                | +              | +   | +                | PA        |
| P36                      | Fresh fruit/cut RTE fruit products       | Melon                                  | -5                | +              | +   | +                | PA        |
| P37                      | Fresh fruit/cut RTE fruit products       | Mango and watermelon fingers           | -5                | +              | +   | +                | PA        |
| P38                      | Fresh fruit/cut RTE fruit products       | Apple pineapple and grape              | -5                | +              | +   | +                | PA        |
| P39                      | Fresh fruit/cut RTE fruit products       | Mango pineapple and passionfruit       | -5                | -              | -   | -                | NA        |
| P40                      | Fresh fruit/cut RTE fruit products       | Melon medley                           | -5                | +              | +   | +                | PA        |





| Fresh produce and fruits |                                   |   |                   |                |   |                  |           |
|--------------------------|-----------------------------------|---|-------------------|----------------|---|------------------|-----------|
| Sample code              | Type                              | Sample used   | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| T41                      | Heat treated fruit and vegetables | Apple and Raspberry juice                               | -1                | -              | -   | -                | NA        |
| T42                      | Heat treated fruit and vegetables | Pure orange juice                                       | -1                | -              | -   | +                | ND        |
| T43                      | Heat treated fruit and vegetables | Apple and pear juice                                    | -1                | -              | -   | -                | NA        |
| T44                      | Heat treated fruit and vegetables | Multivitamin boost fruit juice                          | -1                | -              | -   | -                | NA        |
| T45                      | Heat treated fruit and vegetables | Pure orange and mango fruit juice                       | -1                | +              | +   | +                | PA        |
| P46                      | Heat treated fruit and vegetables | Organic Mangoes Baby Food Pouch                         | -1                | -              | -   | +                | ND        |
| P47                      | Heat treated fruit and vegetables | Organic Squeezy Fruit Pouch Apple Mango & Banana        | -1                | -              | -   | +                | ND        |
| P48                      | Heat treated fruit and vegetables | Organic Sweet Potatoes Food Pouch                       | -1                | +              | +   | +                | PA        |
| P49                      | Heat treated fruit and vegetables | Organic Pears Food Pouch                                | -1                | -              | -   | +                | ND        |
| P50                      | Heat treated fruit and vegetables | Organic Squeezy Fruit Apple Strawberry & Blueberry Food | -1                | -              | -   | -                | NA        |
| P51                      | Heat treated fruit and vegetables | Organic Prunes Baby Food Pouch                          | -1                | -              | -   | -                | NA        |
| P52                      | Heat treated fruit and vegetables | Organic Bananas Food Pouch                              | -1                | -              | -   | -                | NA        |
| P53                      | Heat treated fruit and vegetables | Squeezy Fruit Pouch Apple & Strawberry                  | -1                | -              | -   | -                | NA        |
| P54                      | Heat treated fruit and vegetables | Organic Carrots Baby Food Pouch                         | -1                | -              | -   | -                | NA        |
| P55                      | Heat treated fruit and vegetables | Strawberry & Banana Fruit Smoothie                      | -1                | -              | -   | -                | NA        |
| P56                      | Heat treated fruit and vegetables | Tangy & Sweet Orange Juice with Bits                    | -1                | -              | -   | -                | NA        |
| P57                      | Heat treated fruit and vegetables | Tropical Juice  | -1                | +              | +   | -                | PD        |
| P58                      | Heat treated fruit and vegetables | Pineapple Coconut & Banana Fruit Smoothie               | -1                | -              | -   | -                | NA        |
| P59                      | Heat treated fruit and vegetables | Mango & Passion Fruit Smoothie                          | -1                | -              | -   | -                | NA        |
| P60                      | Heat treated fruit and vegetables | Vegetable juice   | -1                | -              | -   | +                | ND        |



| Multicomponent foods |  |                                      |                   |                |   |                  |           |
|----------------------|--|--------------------------------------|-------------------|----------------|---|------------------|-----------|
| Sample code          | Type   | Sample used                          | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| T46                  | Composite foods with substantial raw ingredients | Chicken triple sandwich              | -7                | +              | +   | +                | PA        |
| T47                  | Composite foods with substantial raw ingredients | Chicken and sweetcorn sandwich       | -7                | +              | +   | +                | PA        |
| T48                  | Composite foods with substantial raw ingredients | Tuna and cucumber sandwich           | -7                | +              | +   | +                | PA        |
| T49                  | Composite foods with substantial raw ingredients | Egg and cress sandwich               | -7                | -              | -   | -                | NA        |
| T50                  | Composite foods with substantial raw ingredients | Cheese and onion sandwich            | -7                | -              | -   | +                | ND        |
| C6                   | Composite foods with substantial raw ingredients | Spicy chicken pasta                  | -7                | -              | -   | -                | NA        |
| C7                   | Composite foods with substantial raw ingredients | Coronation chicken pasta salad       | -7                | -              | -   | -                | NA        |
| C8                   | Composite foods with substantial raw ingredients | Southern fried chicken pasta         | -7                | +              | +   | +                | PA        |
| C9                   | Composite foods with substantial raw ingredients | Onion bhaji sandwich                 | -7                | -              | -   | -                | NA        |
| C10                  | Composite foods with substantial raw ingredients | Spinach pinenut pasta                | -7                | -              | -   | -                | NA        |
| C11                  | Composite foods with substantial raw ingredients | Chicken and bacon pasta              | -7                | -              | -   | -                | NA        |
| C12                  | Composite foods with substantial raw ingredients | Feta and slow roast tomato pasta     | -7                | -              | -   | -                | NA        |
| C13                  | Composite foods with substantial raw ingredients | Cheese and pickle sandwich           | -7                | +              | +   | -                | PD        |
| C14                  | Composite foods with substantial raw ingredients | Chicken and chorizo pasta salad      | -7                | +              | +   | -                | PD        |
| C15                  | Composite foods with substantial raw ingredients | Chicken and bacon pasta              | -7                | +              | +   | +                | PA        |
| C16                  | Composite foods with substantial raw ingredients | Chicken tomato and basil pasta salad | -7                | -              | -   | -                | NA        |
| C17                  | Composite foods with substantial raw ingredients | Cheese and tomato pasta              | -7                | -              | -   | -                | NA        |
| C18                  | Composite foods with substantial raw ingredients | Tomato and basil pasta salad         | -7                | -              | -   | +                | ND        |
| C19                  | Composite foods with substantial raw ingredients | Honey and mustard chicken pasta      | -7                | +              | +   | +                | ND        |



| Multicomponent foods |  |                                  |                   |                |   |                  |           |
|----------------------|--|----------------------------------|-------------------|----------------|---|------------------|-----------|
| Sample code          | Type   | Sample used                      | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| C20                  | Composite foods with substantial raw ingredients | Tuna & sweetcorn pasta           | -7                | +              | +   | +                | PA        |
| T51                  | RTRH/RTE foods (chilled, frozen)                 | Lasagne                          | -6                | +              | +   | +                | PA        |
| T52                  | RTRH/RTE foods (chilled, frozen)                 | Tomato and mozzarella pasta bake | -6                | +              | +   | +                | PA        |
| T53                  | RTRH/RTE foods (chilled, frozen)                 | Macaroni cheese                  | -6                | +              | +   | +                | PA        |
| T54                  | RTRH/RTE foods (chilled, frozen)                 | Tikka masala & pilau rice        | -6                | +              | +   | +                | PA        |
| T55                  | RTRH/RTE foods (chilled, frozen)                 | Spicy chicken pasta              | -6                | -              | -   | -                | NA        |
| C26                  | RTRH/RTE foods (chilled, frozen)                 | Ham & mushroom tagliatelle       | -6                | +              | +   | +                | PA        |
| C27                  | RTRH/RTE foods (chilled, frozen)                 | Chicken hotpot                   | -6                | +              | +   | +                | PA        |
| C28                  | RTRH/RTE foods (chilled, frozen)                 | Cottage pie                      | -6                | +              | +   | +                | PA        |
| C29                  | RTRH/RTE foods (chilled, frozen)                 | Shepards pie                     | -6                | -              | -   | -                | NA        |
| C30                  | RTRH/RTE foods (chilled, frozen)                 | Tomato & mozzarella penne bake   | -6                | +              | +   | +                | PA        |
| C31                  | RTRH/RTE foods (chilled, frozen)                 | Chicken curry & rise             | -6                | -              | -   | -                | NA        |
| C32                  | RTRH/RTE foods (chilled, frozen)                 | Sausage & mash                   | -6                | -              | -   | -                | NA        |
| C33                  | RTRH/RTE foods (chilled, frozen)                 | Spaghetti Bolognese              | -6                | +              | +   | +                | PA        |
| C34                  | RTRH/RTE foods (chilled, frozen)                 | Italian beef lasagne             | -6                | -              | -   | -                | NA        |
| C35                  | RTRH/RTE foods (chilled, frozen)                 | Chicken korma with pilau rice    | -6                | +              | +   | +                | PA        |
| C36                  | RTRH/RTE foods (chilled, frozen)                 | Beef stew & dumplings            | -6                | +              | +   | +                | PA        |
| C37                  | RTRH/RTE foods (chilled, frozen)                 | Chicken & bacon pasta bake       | -6                | +              | +   | +                | PA        |
| C38                  | RTRH/RTE foods (chilled, frozen)                 | Spaghetti & meatballs            | -6                | -              | -   | -                | NA        |
| C39                  | RTRH/RTE foods (chilled, frozen)                 | Spaghetti carbonara              | -6                | +              | +   | +                | PA        |



| Multicomponent foods |                                  |  |                   |                |   |                  |           |
|----------------------|----------------------------------|--|-------------------|----------------|---|------------------|-----------|
| Sample code          | Type                             | Sample used                            | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| C40                  | RTRH/RTE foods (chilled, frozen) | Macaroni and cheese                    | -6                | -              | -   | -                | NA        |
| T56                  | Mayonnaise based deli salads     | Triple grain salad mayonnaise          | -8                | -              | -   | -                | NA        |
| T57                  | Mayonnaise based deli salads     | Party salad mayonnaise                 | -8                | -              | -   | -                | NA        |
| T58                  | Mayonnaise based deli salads     | Chicken Caesar salad                   | -8                | -              | -   | +                | ND        |
| T59                  | Mayonnaise based deli salads     | Ham egg and coleslaw salad             | -8                | -              | -   | -                | NA        |
| T60                  | Mayonnaise based deli salads     | Chicken & bacon salad                  | -8                | +              | +   | +                | PA        |
| C45                  | Mayonnaise based deli salads     | Cesear salad                           | -8                | +              | +   | +                | PA        |
| C47                  | Mayonnaise based deli salads     | Chicken and bacon Caesar salad         | -8                | +              | +   | +                | PA        |
| C48                  | Mayonnaise based deli salads     | Celery fruit nut salad                 | -8                | -              | -   | -                | NA        |
| C49                  | Mayonnaise based deli salads     | Tuna and sweetcorn salad               | -8                | -              | -   | -                | NA        |
| C50                  | Mayonnaise based deli salads     | Falafel and rice salad                 | -8                | -              | -   | -                | NA        |
| C51                  | Mayonnaise based deli salads     | Chicken rice bowl                      | -8                | -              | -   | -                | NA        |
| C52                  | Mayonnaise based deli salads     | Bang bang chicken salad                | -8                | +              | +   | +                | PA        |
| C53                  | Mayonnaise based deli salads     | Chicken, broccoli and almond salad     | -8                | +              | +   | +                | PA        |
| C54                  | Mayonnaise based deli salads     | Hummus salad                           | -8                | -              | -   | -                | NA        |
| C55                  | Mayonnaise based deli salads     | Prawn layered salad                    | -8                | -              | -   | -                | NA        |
| C56                  | Mayonnaise based deli salads     | Baby potato & free-range egg salad     | -8                | -              | -   | -                | NA        |
| C57                  | Mayonnaise based deli salads     | Chicken & bacon Caesar pasta salad     | -8                | -              | -   | -                | NA        |
| C58                  | Mayonnaise based deli salads     | Baby potato salad                      | -8                | -              | -   | -                | NA        |
| C59                  | Mayonnaise based deli salads     | Cauliflower & spicy rice w/ mayonnaise | -8                | +              | +   | +                | PA        |

2022LR113 Neogen Soleris Non-Fermenting Total  
Viable Count (NF-TVC) vial MCS ILS Summary  
Report 14/10/2024



| Multicomponent foods |                              |                             |                   |                |   |                  |           |
|----------------------|------------------------------|-----------------------------|-------------------|----------------|---|------------------|-----------|
| Sample code          | Type                         | Sample used                 | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| C60                  | Mayonnaise based deli salads | Chicken & bacon pasta salad | -8                | -              | -   | -                | NA        |



| Raw and Ready to Cook Meat and Poultry |                           |   |                   |                |   |                  |           |
|--|---------------------------|---|-------------------|----------------|---|------------------|-----------|
| Sample code                            | Type                      | Sample used                                     | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| T1                                     | Raw poultry and meat cuts | Chicken breast fillets                          | -7                | -              | -   | -                | NA        |
| T2                                     | Raw poultry and meat cuts | Lean diced beef                                 | -7                | -              | -   | -                | NA        |
| T3                                     | Raw poultry and meat cuts | Fresh lamb chops                                | -7                | +              | +   | +                | PA        |
| T4                                     | Raw poultry and meat cuts | Turkey thigh mince 7% fat                       | -7                | +              | +   | -                | PD        |
| T5                                     | Raw poultry and meat cuts | Pork loin steaks                                | -7                | -              | -   | -                | NA        |
| R6                                     | Raw poultry and meat cuts | Turkey breast fillet portions                   | -7                | +              | +   | +                | PA        |
| R7                                     | Raw poultry and meat cuts | Turkey burgers                                  | -7                | +              | +   | +                | PA        |
| R8                                     | Raw poultry and meat cuts | Turkey breast steaks                            | -7                | -              | -   | -                | NA        |
| R9                                     | Raw poultry and meat cuts | Turkey streaks                                  | -7                | -              | -   | -                | NA        |
| R10                                    | Raw poultry and meat cuts | Turkey breast mini fillets                      | -7                | -              | -   | +                | ND        |
| R11                                    | Raw poultry and meat cuts | Chicken skin on drumsticks                      | -7                | +              | +   | +                | PA        |
| R12                                    | Raw poultry and meat cuts | Fresh small whole chicken                       | -7                | -              | -   | -                | NA        |
| R13                                    | Raw poultry and meat cuts | Fresh chicken thigh fillets skinless + boneless | -7                | -              | -   | -                | NA        |
| R14                                    | Raw poultry and meat cuts | Chicken breast mini fillets                     | -7                | -              | -   | -                | NA        |
| R15                                    | Raw poultry and meat cuts | Chicken skin on thighs                          | -7                | +              | +   | +                | PA        |
| R16                                    | Raw poultry and meat cuts | Chicken breast 2                                | -7                | -              | -   | -                | NA        |
| R17                                    | Raw poultry and meat cuts | Chicken breast fillets 3                        | -7                | +              | +   | +                | PA        |
| R18                                    | Raw poultry and meat cuts | British pork shoulder steaks                    | -7                | -              | -   | -                | NA        |
| R19                                    | Raw poultry and meat cuts | Thick cut british pork chops                    | -7                | -              | -   | -                | NA        |
| R20                                    | Raw poultry and meat cuts | Diced beef                                      | -7                | +              | +   | +                | PA        |



| Raw and Ready to Cook Meat and Poultry |                    |                                |                   |                |   |                  |           |
|--|--------------------|--------------------------------|-------------------|----------------|---|------------------|-----------|
| Sample code                            | Type               | Sample used                    | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| T6                                     | Raw processed meat | BBQ pork riblets               | -7                | -              | -   | -                | NA        |
| T7                                     | Raw processed meat | Pork shoulder in a bbq sauce   | -7                | -              | -   | +                | ND        |
| T8                                     | Raw processed meat | Sweet and smoky beef kebabs    | -7                | +              | +   | -                | PD        |
| T9                                     | Raw processed meat | Beef burgers                   | -7                | +              | +   | +                | PA        |
| T10                                    | Raw processed meat | Pork sausages                  | -7                | +              | +   | +                | PA        |
| R26                                    | Raw processed meat | Turkey gluten free sausages    | -7                | -              | -   | -                | NA        |
| R27                                    | Raw processed meat | 2% fat turkey mince            | -7                | +              | +   | +                | PA        |
| R28                                    | Raw processed meat | Turkey meatballs               | -7                | +              | +   | +                | PA        |
| R29                                    | Raw processed meat | British pork stir fry          | -7                | -              | -   | -                | NA        |
| R30                                    | Raw processed meat | British pork belly slices      | -7                | -              | -   | -                | NA        |
| R31                                    | Raw processed meat | Pork tenderloin fillet         | -7                | +              | +   | +                | PA        |
| R32                                    | Raw processed meat | Smoked gammon steaks           | -7                | +              | +   | +                | PA        |
| R33                                    | Raw processed meat | Slow cooked pulled pork        | -7                | -              | -   | -                | NA        |
| R34                                    | Raw processed meat | Chinese style pork             | -7                | -              | -   | -                | NA        |
| R35                                    | Raw processed meat | 5% fat beef mince              | -7                | +              | +   | +                | PA        |
| R36                                    | Raw processed meat | 10% fat meatballs              | -7                | -              | -   | -                | NA        |
| R37                                    | Raw processed meat | Quarter pounder burgers        | -7                | +              | +   | -                | PD        |
| R38                                    | Raw processed meat | Caramelised onion beef burgers | -7                | -              | -   | -                | NA        |
| R39                                    | Raw processed meat | Wagyu beef burgers             | -7                | -              | -   | -                | NA        |
| R40                                    | Raw processed meat | Venison burgers                | -7                | +              | +   | -                | PD        |



| Raw and Ready to Cook Meat and Poultry |                      |   |                   |                |   |                  |           |
|--|----------------------|---|-------------------|----------------|---|------------------|-----------|
| Sample code                            | Type                 | Sample used                                 | Dilution selected | Soleris result | PCA confirmation result from Soleris vial | Reference result | Agreement |
| T11                                    | RTC meat and poultry | Breaded chicken goujons                     | -4                | -              | -   | -                | NA        |
| T12                                    | RTC meat and poultry | Turkey meatballs                            | -4                | -              | -   | -                | NA        |
| T13                                    | RTC meat and poultry | BBQ roast chicken wings                     | -4                | +              | +   | +                | PA        |
| T14                                    | RTC meat and poultry | Chicken kiev bites                          | -4                | +              | +   | +                | PA        |
| T15                                    | RTC meat and poultry | Southern fried breaded chicken mini fillets | -4                | -              | -   | -                | NA        |
| R46                                    | RTC meat and poultry | Chicken quarter pounder                     | -4                | -              | -   | -                | NA        |
| R47                                    | RTC meat and poultry | Pork & herb sausage patties                 | -4                | +              | +   | +                | PA        |
| R48                                    | RTC meat and poultry | Garlic chicken kiev                         | -4                | +              | +   | +                | PA        |
| R49                                    | RTC meat and poultry | Garlic butter chicken kiev                  | -4                | +              | +   | +                | PA        |
| R50                                    | RTC meat and poultry | Buttered chicken breast joint               | -4                | +              | +   | +                | PA        |
| R51                                    | RTC meat and poultry | Chicken wings side                          | -4                | +              | +   | +                | PA        |
| R52                                    | RTC meat and poultry | Wild garlic chicken kiev with salted butter | -4                | +              | +   | +                | PA        |
| R53                                    | RTC meat and poultry | Breaded garlic chicken kiev                 | -4                | +              | +   | +                | PA        |
| R54                                    | RTC meat and poultry | Garlic breast kiev whole fillet             | -4                | +              | +   | +                | PA        |
| R55                                    | RTC meat and poultry | Cheese & bacon chicken kiev                 | -4                | +              | +   | +                | PA        |
| R56                                    | RTC meat and poultry | Beef & herb meatballs                       | -4                | -              | -   | -                | NA        |
| R57                                    | RTC meat and poultry | Beef meatballs                              | -4                | -              | -   | -                | NA        |
| R58                                    | RTC meat and poultry | British beef quarterpounders                | -4                | +              | +   | -                | PD        |
| R59                                    | RTC meat and poultry | Turkey mini kiev                            | -4                | +              | +   | +                | PA        |
| R60                                    | RTC meat and poultry | Cheese and bacon chicken kiev               | -4                | +              | +   | +                | PA        |



**ANNEX D: Raw data from the RLOD study**

| Sample code                                | Level | Inoculation (cfu/g) | Soleris alternative method      |                           |                                | Reference method |
|--|-------|---------------------|---------------------------------|---------------------------|--------------------------------|------------------|
|  |       |                     | Soleris result (detection time) | Confirmation result (+/-) | Alternative final result (+/-) |                  |
| Category - Dairy products (heat processed) |       |                     |                                 |                           |                                |                  |
| R1   | Blank | N/A                 | ND                              | -                         | -                              | -                |
| R2   |       |                     | ND                              | -                         | -                              | -                |
| R3   |       |                     | ND                              | -                         | -                              | -                |
| R4   |       |                     | ND                              | -                         | -                              | -                |
| R5   |       |                     | ND                              | -                         | -                              | -                |
| R6   | Low   | 11                  | ND                              | -                         | -                              | +                |
| R7   |       |                     | ND                              | -                         | -                              | -                |
| R8   |       |                     | ND                              | +                         | -                              | +                |
| R9   |       |                     | 10.3                            | +                         | +                              | +                |
| R10  |       |                     | 14.3                            | +                         | +                              | +                |
| R11  |       |                     | 11.6                            | +                         | +                              | -                |
| R12  |       |                     | 11.7                            | +                         | +                              | -                |
| R13  |       |                     | ND                              | -                         | -                              | +                |
| R14  |       |                     | ND                              | -                         | -                              | -                |
| R15  |       |                     | 10.4                            | +                         | +                              | +                |
| R16  |       |                     | 11.4                            | +                         | +                              | +                |
| R17  |       |                     | 13                              | +                         | +                              | -                |
| R18  |       |                     | 10.4                            | +                         | +                              | -                |
| R19  |       |                     | 14.3                            | +                         | +                              | -                |
| R20  |       |                     | ND                              | +                         | -                              | +                |
| R21  |       |                     | ND                              | -                         | -                              | +                |
| R22  |       |                     | 11.2                            | +                         | +                              | +                |
| R23  |       |                     | 10.1                            | +                         | +                              | +                |
| R24  |       |                     | 13.5                            | +                         | +                              | +                |
| R25  |       |                     | 11.7                            | +                         | +                              | +                |
| R26  | High  | 25                  | 12.7                            | +                         | +                              | +                |
| R27  |       |                     | 11.7                            | +                         | +                              | -                |
| R28  |       |                     | 23.7                            | +                         | +                              | +                |
| R29  |       |                     | 10.9                            | +                         | +                              | +                |
| R30  |       |                     | 12.6                            | +                         | +                              | +                |

| Sample code                             | Level | Inoculation (cfu/g) | Soleris alternative method      |                           |                                | Reference method |
|---|-------|---------------------|---------------------------------|---------------------------|--------------------------------|------------------|
|   |       |                     | Soleris result (detection time) | Confirmation result (+/-) | Alternative final result (+/-) |                  |
| Category - Raw and RTC fishery products |       |                     |                                 |                           |                                |                  |
| R31                                     | Blank | N/A                 | ND                              | -                         | -                              | -                |
| R32                                     |       |                     | ND                              | -                         | -                              | -                |
| R33                                     |       |                     | ND                              | -                         | -                              | -                |
| R34                                     |       |                     | ND                              | -                         | -                              | -                |
| R35                                     |       |                     | ND                              | -                         | -                              | -                |
| R36                                     | Low   | 12                  | ND                              | -                         | -                              | -                |
| R37                                     |       |                     | ND                              | -                         | -                              | -                |
| R38                                     |       |                     | ND                              | -                         | -                              | -                |
| R39                                     |       |                     | ND                              | -                         | -                              | -                |
| R40                                     |       |                     | 25.7                            | -                         | -                              | +                |
| R41                                     |       |                     | 13.4                            | +                         | +                              | +                |
| R42                                     |       |                     | 22.8                            | +                         | +                              | +                |
| R43                                     |       |                     | 13.5                            | +                         | +                              | +                |
| R44                                     |       |                     | 14                              | +                         | +                              | +                |
| R45                                     |       |                     | 14.9                            | +                         | +                              | +                |
| R46                                     |       |                     | 13.9                            | +                         | +                              | +                |
| R47                                     |       |                     | 13.2                            | +                         | +                              | +                |
| R48                                     |       |                     | 14.7                            | +                         | +                              | +                |
| R49                                     |       |                     | 26.3                            | -                         | -                              | +                |
| R50                                     |       |                     | 20.8                            | -                         | -                              | +                |
| R51                                     |       |                     | 12.3                            | +                         | +                              | +                |
| R52                                     |       |                     | ND                              | -                         | -                              | +                |
| R53                                     |       |                     | 13.5                            | +                         | +                              | -                |
| R54                                     |       |                     | 15.1                            | +                         | +                              | +                |
| R55                                     |       |                     | 14.6                            | +                         | +                              | +                |
| R56                                     | High  | 24                  | 12.8                            | +                         | +                              | +                |
| R57                                     |       |                     | 13.1                            | +                         | +                              | +                |
| R58                                     |       |                     | 13.2                            | +                         | +                              | +                |
| R59                                     |       |                     | 13.1                            | +                         | +                              | +                |
| R60                                     |       |                     | 13.1                            | +                         | +                              | +                |

| Sample code  | Level | Inoculation (cfu/g) | Soleris alternative method      |                           |                                | Reference method |
|--|-------|---------------------|---------------------------------|---------------------------|--------------------------------|------------------|
|  |       |                     | Soleris result (detection time) | Confirmation result (+/-) | Alternative final result (+/-) |                  |
| Category - Produce and fruits (combined category: fresh and processed) |       |                     |                                 |                           |                                |                  |
| R61  | Blank | N/A                 | ND                              | -                         | -                              | -                |
| R62  |       |                     | ND                              | -                         | -                              | -                |
| R63  |       |                     | ND                              | -                         | -                              | -                |
| R64  |       |                     | ND                              | -                         | -                              | -                |
| R65  |       |                     | ND                              | -                         | -                              | -                |
| R66  | Low   | 11                  | 15.4                            | +                         | +                              | +                |
| R67  |       |                     | 15                              | +                         | +                              | +                |
| R68  |       |                     | 15.7                            | +                         | +                              | +                |
| R69  |       |                     | 17                              | +                         | +                              | +                |
| R70  |       |                     | 15.4                            | +                         | +                              | +                |
| R71  |       |                     | 16.9                            | +                         | +                              | +                |
| R72  |       |                     | 15.1                            | +                         | +                              | +                |
| R73  |       |                     | 15.4                            | +                         | +                              | +                |
| R74  |       |                     | 15.6                            | +                         | +                              | +                |
| R75  |       |                     | 15.7                            | +                         | +                              | +                |
| R76  |       |                     | 14.8                            | +                         | +                              | +                |
| R77  |       |                     | 16.1                            | +                         | +                              | +                |
| R78  |       |                     | 15.9                            | +                         | +                              | +                |
| R79  |       |                     | ND                              |                           | -                              | +                |
| R80  |       |                     | 16.4                            | +                         | +                              | -                |
| R81  |       |                     | ND                              | -                         | -                              | -                |
| R82  |       |                     | 15.5                            | +                         | +                              | -                |
| R83  |       |                     | ND                              | -                         | -                              | +                |
| R84  |       |                     | ND                              | -                         | -                              | -                |
| R85  |       |                     | ND                              | -                         | -                              | +                |
| R86  | High  | 23                  | 15                              | +                         | +                              | +                |
| R87  |       |                     | 15.9                            | +                         | +                              | +                |
| R88  |       |                     | 15.5                            | +                         | +                              | +                |
| R89  |       |                     | 17.3                            | +                         | +                              | +                |
| R90  |       |                     | 16.4                            | +                         | +                              | +                |

| Sample code  | Level | Inoculation (cfu/g) | Soleris alternative method      |                           |                                | Reference method |
|--|-------|---------------------|---------------------------------|---------------------------|--------------------------------|------------------|
|  |       |                     | Soleris result (detection time) | Confirmation result (+/-) | Alternative final result (+/-) |                  |
| Category - Raw and RTC meat and poultry (combined) |       |                     |                                 |                           |                                |                  |
| R91  | Blank | N/A                 | ND                              | -                         | -                              | -                |
| R92  |       |                     | ND                              | -                         | -                              | -                |
| R93  |       |                     | ND                              | -                         | -                              | -                |
| R94  |       |                     | ND                              | -                         | -                              | -                |
| R95  |       |                     | ND                              | -                         | -                              | -                |
| R96  | Low   | 11                  | ND                              | -                         | -                              | +                |
| R97  |       |                     | ND                              | -                         | -                              | -                |
| R98  |       |                     | ND                              | -                         | -                              | -                |
| R99  |       |                     | ND                              | -                         | -                              | -                |
| R100   |       |                     | ND                              | -                         | -                              | -                |
| R101   |       |                     | 23.4                            | +                         | +                              | -                |
| R102   |       |                     | 23.3                            | +                         | +                              | +                |
| R103   |       |                     | 21.7                            | +                         | +                              | -                |
| R104   |       |                     | 26.3                            | +                         | +                              | -                |
| R105   |       |                     | 25.6                            | +                         | +                              | +                |
| R106   |       |                     | 22.8                            | +                         | +                              | +                |
| R107   |       |                     | ND                              | -                         | -                              | +                |
| R108   |       |                     | ND                              | -                         | -                              | -                |
| R109   |       |                     | 25.2                            | +                         | +                              | -                |
| R110   |       |                     | ND                              | -                         | -                              | +                |
| R111   |       |                     | 24.1                            | +                         | +                              | -                |
| R112   |       |                     | 29.5                            | +                         | +                              | -                |
| R113   |       |                     | 25.6                            | +                         | +                              | +                |
| R114   |       |                     | 22.5                            | +                         | +                              | +                |
| R115   |       |                     | ND                              | -                         | -                              | -                |
| R116   | High  | 22                  | 23.2                            | +                         | +                              | +                |
| R117   |       |                     | 22.4                            | +                         | +                              | +                |
| R118   |       |                     | 23                              | +                         | +                              | +                |
| R119   |       |                     | 22.5                            | +                         | +                              | +                |
| R120   |       |                     | 22.6                            | +                         | +                              | +                |

| Sample code  | Level | Inoculation (cfu/g) | Soleris alternative method      |                           |                                | Reference method |
|--|-------|---------------------|---------------------------------|---------------------------|--------------------------------|------------------|
|  |       |                     | Soleris result (detection time) | Confirmation result (+/-) | Alternative final result (+/-) |                  |
| Category - Multicomponent foods or meal components |       |                     |                                 |                           |                                |                  |
| R121   | Blank | N/A                 | ND                              | -                         | -                              | -                |
| R122   |       |                     | ND                              | -                         | -                              | -                |
| R123   |       |                     | ND                              | -                         | -                              | -                |
| R124   |       |                     | ND                              | -                         | -                              | -                |
| R125   |       |                     | ND                              | -                         | -                              | -                |
| R126   | Low   | 11                  | 17.5                            | +                         | +                              | +                |
| R127   |       |                     | 18.4                            | +                         | +                              | +                |
| R128   |       |                     | 16.8                            | +                         | +                              | -                |
| R129   |       |                     | 19.4                            | +                         | +                              | +                |
| R130   |       |                     | 17.7                            | +                         | +                              | +                |
| R131   |       |                     | 17.5                            | +                         | +                              | +                |
| R132   |       |                     | 15.5                            | +                         | +                              | +                |
| R133   |       |                     | ND                              | -                         | -                              | -                |
| R134   |       |                     | ND                              | -                         | -                              | -                |
| R135   |       |                     | 16.1                            | +                         | +                              | +                |
| R136   |       |                     | 16.9                            | +                         | +                              | +                |
| R137   |       |                     | 17.5                            | +                         | +                              | +                |
| R138   |       |                     | 20.3                            | +                         | +                              | -                |
| R139   |       |                     | 17.1                            | +                         | +                              | +                |
| R140   |       |                     | 15.8                            | +                         | +                              | -                |
| R141   |       |                     | 16                              | +                         | +                              | +                |
| R142   |       |                     | 16.4                            | +                         | +                              | -                |
| R143   |       |                     | 17.1                            | +                         | +                              | +                |
| R144   |       |                     | 17.3                            | +                         | +                              | -                |
| R145   |       |                     | 16.8                            | +                         | +                              | +                |
| R146   | High  | 28                  | 15.9                            | +                         | +                              | +                |
| R147   |       |                     | 15.7                            | +                         | +                              | +                |
| R148   |       |                     | 17.1                            | +                         | +                              | +                |
| R149   |       |                     | 16.1                            | +                         | +                              | +                |
| R150   |       |                     | 16.8                            | +                         | +                              | +                |



# ANNEX E: Raw data on inclusivity and exclusivity

| Code | Genus                   | species                   | CRA Number          | Origin                         | Cfu/plate | Alternative method result |                    |
|------|-------------------------|---------------------------|---------------------|--------------------------------|-----------|---------------------------|--------------------|
|      |                         |                           |                     |                                |           | Soleris result            | Detection time (h) |
| 1    | <i>Raoultella</i>       | <i>terrigena</i>          | 17343               | raw milk                       | 30        | D                         | 14.1               |
| 2    | <i>Enterobacter</i>     | <i>cloacae</i>            | 1472                | dried milk                     | 108       | D                         | 10.9               |
| 3    | <i>Klebsiella</i>       | <i>oxytoca</i>            | 8387                | Water                          | 10        | D                         | 14.7               |
| 4    | <i>Kluyvera</i>         | <i>ascorbata</i>          | 17126               | industrial                     | 20        | D                         | 14.6               |
| 5    | <i>Escherichia</i>      | <i>adecarboxylata</i>     | 5501                | Skim milk powder               | 23        | D                         | 12                 |
| 6    | <i>Klebsiella</i>       | <i>trevisanii</i>         | NCIMB 8606          | Ropy cream                     | 18        | D                         | 12.8               |
| 7    | <i>Pantoea</i>          | <i>agglomerans</i>        | 17030, NCIMB 702072 | Pasteurised milk               | 4         | D                         | 19.7               |
| 8    | <i>Aeromonas</i>        | <i>salmonicida</i>        | 8388, NCTC 8049     | tin of milk with a fishy odour | 2         | D                         | 22.9               |
| 9    | <i>Escherichia</i>      | <i>coli</i>               | 1476                | Dried milk                     | 31        | D                         | 12.6               |
| 10   | <i>Rahnella</i>         | <i>aqualtilis</i>         | 16911               | drinking water                 | 40        | D                         | 19.3               |
| 11   | <i>Bacillus</i>         | <i>weihenstephanensis</i> | 16578               | Pasteurised milk               | 3         | D                         | 15.6               |
| 12   | <i>Lysinibacillus</i>   | <i>sphaericus</i>         | 7746                | unknown                        | 1         | D                         | 15.1               |
| 13   | <i>Bacillus</i>         | <i>coagulans</i>          | 16586               | Sterilised milk                | 51        | D                         | 28.6               |
| 14   | <i>Clostridium</i>      | <i>perfringens</i>        | 15911, NCTC 8239    | salt beef                      | 129       | D                         | 16.4               |
| 15   | <i>Staphylococcus</i>   | <i>delphini</i>           | 16891               | factory isolate                | 11        | D                         | 16.8               |
| 16   | <i>Leuconostoc</i>      | <i>mesenteroides</i>      | 17490               | Green ham                      | 51        | D                         | 12.6               |
| 17   | <i>Buttiauxella</i>     | <i>agrestis</i>           | 17110               | Pond water                     | 72        | D                         | 14.2               |
| 18   | <i>Citrobacter</i>      | <i>youngae</i>            | NCTC 13709          | Meat scraps                    | 87        | D                         | 12.3               |
| 19   | <i>Moraxella</i>        | <i>osloensis</i>          | 17043               | milk                           | 46        | D                         | 24.5               |
| 20   | <i>Paraburkholderia</i> | <i>cepaciae</i>           | 16779               | soft drinks environment        | 68        | D                         | 12.8               |
| 21   | <i>Staphylococcus</i>   | <i>carneus</i>            | 284                 | goat's milk                    | 20        | D                         | 13                 |
| 22   | <i>Listeria</i>         | <i>ivanovii</i>           | 1123                | soft cheese                    | 9         | D                         | 19.8               |
| 23   | <i>Streptococcus</i>    | <i>thermophilus</i>       | 16045, NCIMB 8510   | Pasteurised milk               | 20        | D                         | 15.8               |
| 24   | <i>Lactobacillus</i>    | <i>acidophilus</i>        | 7675                | Dairy product                  | 30        | D                         | 16.5               |
| 25   | <i>Carnobacterium</i>   | <i>divergens</i>          | 3910                | Brie                           | 15        | D                         | 15                 |



| Code | Genus                 | species                                    | CRA Number | Origin                        | Cfu/plate | Alternative method result |                    |
|------|-----------------------|--|------------|-------------------------------|-----------|---------------------------|--------------------|
|      |                       |  |            |                               |           | Soleris result            | Detection time (h) |
| 26   | <i>Staphylococcus</i> | <i>saprophyticus</i>                       | 8999       | distilled water environmental | 34        | D                         | 18.9               |
| 27   | <i>Micrococcus</i>    | <i>luteus</i>                              | 3503       | Air sample                    | 8         | D                         | 16.7               |
| 28   | <i>Enterococcus</i>   | <i>faecalis</i>                            | 1513       | Dried milk powder             | 40        | D                         | 15.1               |
| 29   | <i>Staphylococcus</i> | <i>cohnii</i>                              | 272        | skin                          | 83        | D                         | 25.1               |
| 30   | <i>Enterococcus</i>   | <i>faecium</i>                             | 16866      | Uncooked Sausage              | 65        | D                         | 14.9               |
| 31   | <i>Staphylococcus</i> | <i>aureus</i>                              | 409/3026   | Slow cheese                   | 52        | D                         | 14.9               |
| 32   | <i>Staphylococcus</i> | <i>epidermidis</i>                         | 314        | runway & can seam             | 7         | D                         | 20.9               |
| 33   | <i>Pediococcus</i>    | <i>pentosaceus</i>                         | 16030      | Brine                         | 85        | D                         | 28.5               |
| 34   | <i>Listeria</i>       | <i>monocytogenes 1/2a</i>                  | 1100       | Stilton                       | 50        | D                         | 22.1               |
| 35   | <i>Listeria</i>       | <i>innocua</i>                             | 3130       | Cheese factory                | 39        | D                         | 22                 |
| 36   | <i>Listeria</i>       | <i>fleischmanii subspecies fleishmanii</i> | 16876      | Swiss hard cheese             | 105       | D                         | 21.9               |
| 37   | <i>Streptococcus</i>  | <i>parauberis</i>                          | 17214      | Raw mince                     | 31        | D                         | 11.9               |
| 38   | <i>Staphylococcus</i> | <i>hominis</i>                             | 16828      | unknown                       | 5         | D                         | 21.1               |
| 39   | <i>Staphylococcus</i> | <i>warneri</i>                             | 3198       | Dry sausage                   | 11        | D                         | 22.4               |
| 40   | <i>Lactococcus</i>    | <i>lactis</i>                              | 16029      | Green ham                     | 35        | D                         | 25.4               |
| 41   | <i>Micrococcus</i>    | <i>roseus</i>                              | 7775       | water                         | 14        | D                         | 25.5               |
| 42   | <i>Streptococcus</i>  | <i>lactis</i>                              | 1511       | dried milk powder             | 7         | D                         | 16.6               |
| 43   | <i>Enterococcus</i>   | <i>pseudoavium</i>                         | 16852      | Cow udder - bovine mastitis   | 25        | D                         | 13.9               |
| 44   | <i>Candida</i>        | <i>krussei</i>                             | CRA629     | Yogurt base                   | 1         | D                         | 26.6               |
| 45   | <i>Kluyveromyces</i>  | <i>marxianus</i>                           | CRA 6749   | Dairy isolate                 | 25        | D                         | 16.4               |
| 46   | <i>Torulaspora</i>    | <i>delbrukeii</i>                          | CRA16154   | Spoiled yogurt                | 13        | D                         | 24.2               |
| 47   | <i>Geochium</i>       | <i>condidum</i>                            | 14398      | Factory isolate               | 7         | D                         | 16.2               |
| 48   | <i>Fusarium</i>       | <i>solani</i>                              | 16976      | Factory isolate               | 25        | D                         | 26.9               |
| 49   | <i>Debracmyces</i>    | <i>hansenii</i>                            | 15969      | Environment                   | 49        | D                         | 29.6               |
| 50   | <i>Kluyveromyces</i>  | <i>lactis</i>                              | MUCL 28769 | gassy cheese                  | 12        | D                         | 14.4               |

Key D = detected

## ANNEX F: Collaborators in ILS

| Laboratories                                   | Country       | Address                 | Number of collaborators |
|--|---------------|-------------------------|-------------------------|
| NQAC Konolfingen                               | Switzerland   | Bern, Switzerland       | 1                       |
| Aliminter                                      | Spain         | Murcia, Spain           | 1                       |
| Kraft Heinz Elst                               | Netherlands   | Elst, Netherlands       | 2                       |
| Dairygold                                      | Ireland       | Cork, Ireland           | 1                       |
| Danone - Aqualab                               | France        | Évian, France           | 2                       |
| Q labs   | United States | Ohio, US                | 1                       |
| Muller Yoghurts & Desserts                     | England       | Market Drayton, England | 2                       |
| Kraft Heinz Kitt Green                         | England       | Wigan, England          | 2                       |
| PZ Cussons                                     | England       | Manchester, England     | 2                       |
| Campden BRI – Microbiology Analytical Services | England       | Gloucestershire, UK     | 1                       |